

Cryogenic Vials with Green Caps, 1.2 mL



Sterile, self-standing, polypropylene cryogenic vials with green caps

Scientists Helping Scientists™ | WWW.STEMCELL.COM

Catalog # 100-0095 125 Vials
38047 500 Vials

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713
INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM
FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Cryogenic Vials with Green Caps are ideal for reliable sample cryostorage. The sterilized, clear polypropylene vial withstands temperatures as low as -196°C (in gas/vapor phase) and has a self-standing design. The green polypropylene cap includes a silicone washer for a dependable seal. The vial features graduated markings and a frosted side panel for easy labeling.

Catalog #100-0095, 1.2 mL - 125 Vials; 125 Vials/pack (Trueline Catalog #TR7000 and equivalent to Corning Catalog #430487)
Catalog #38047, 1.2 mL - 500 Vials; 125 Vials/pack (Trueline Catalog #TR7000 and equivalent to Corning Catalog #430487)

APPLICATIONS

- Cryopreservation
- Cell storage

Properties

Storage: Store at 15 - 25°C.

Shelf Life: For Catalog #100-0095: Stable until expiry date (EXP) on STEMCELL label.
For Catalog #38047: Stable until expiry date on manufacturer label.

Features

MATERIALS

- Tube: clear polypropylene, conical-bottom, self-standing
- Cap: green polypropylene, internally threaded
- Washer: silicone
- Withstands temperatures as low as -196°C (in gas/vapor phase)
- RNase/DNase-free and non-pyrogenic

DIMENSIONS

- Length - Vial with cap: 41 mm
- Diameter - Vial: 12.7 mm; Cap: 11.9 mm

NOMINAL VOLUME

- 1.2 mL
- Graduations at 0.5 mL and 1 mL

STERILITY

- Gamma irradiation

Handling / Directions For Use

Use as directed in protocol of choice.

Notes and Tips

For chemical resistance or thermoplastic property information, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products available from STEMCELL Technologies, visit www.stemcell.com or contact us at techsupport@stemcell.com.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2019 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.