

35 mm Dish, Tissue Culture-Treated



Sterile 35 mm tissue culture-treated dish

Catalog # 100-0101 20 Dishes
38050 500 Dishes

Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Tissue Culture-Treated Dishes are ideal for attachment-dependent cell and tissue culture. These polystyrene dishes are transparent for superior optical clarity and are stackable for easy handling and storage. The dishes have been sterilized and are vented to allow for gas exchange.

Catalog #100-0101, 35 mm - 20 Dishes; 20 Dishes/bag (Trueline Catalog #TR4000 and equivalent to Corning Catalog #353001/430165)
Catalog #38050, 35 mm - 500 Dishes; 20 Dishes/bag (Trueline Catalog #TR4000 and equivalent to Corning Catalog #353001/430165)

APPLICATIONS

- Cell culture
- Microscopy

Properties

Storage: Store at 15 - 25°C.

Shelf Life: For Catalog #100-0101: Stable until expiry date (EXP) on STEMCELL label.

For Catalog #38050: Stable until expiry date on manufacturer label.

Features

MATERIALS

- Dish and lid: clear and transparent polystyrene
- Tissue culture-treated for optimal cell attachment
- Non-pyrogenic

DIMENSIONS

- 35 x 10 mm (diameter x height)
- Stackable for easy handling and storage

CAPACITY

- Recommended working volume of 1.8 - 2.7 mL
- Approximate cell growth area of 8 cm²

STERILITY

- Gamma irradiation

Handling / Directions For Use

Use as directed in protocol of choice.

Notes and Tips

For chemical resistance or thermoplastic property information, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products available from STEMCELL Technologies, visit www.stemcell.com or contact us at techsupport@stemcell.com.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2019 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.