

Round-Bottom Polypropylene Tubes without Caps, 5 mL

Polypropylene tubes without caps

Catalog # 100-0094
38056

125 Tubes
1000 Tubes



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Falcon® Round-Bottom Polypropylene Tubes, 5 mL are sterile, high-quality tubes ideal for a wide variety of sample preparation, sample storage, and sample centrifugation tasks. They are required for select protocols in cell isolation, flow cytometry, and fluorescence-activated cell sorting (FACS). Tubes are RNase-/DNase-free.

Catalog #100-0094, 5 mL - 125 Tubes without caps; 125 Tubes/pack (Corning Catalog #352053)

Catalog #38056, 5mL - 1000 Tubes without caps; 125 Tubes/pack (Corning Catalog #352053)

APPLICATIONS

- Cell separation
- Cell culture
- Sample storage
- Liquid handling

Properties

Storage: Store at 15 - 25°C.

Shelf Life: For Catalog #100-0094: Stable until expiry date (EXP) on STEMCELL label.

For Catalog #38056: Stable until expiry date on manufacturer label.

Features

MATERIALS

- Clear polypropylene
- Biologically inert, RNase-/DNase-free, and non-pyrogenic

DIMENSIONS

- 12 x 74 mm (outer diameter x height)

NOMINAL VOLUME

- 5 mL
- No graduations

CENTRIFUGAL FORCE

- Up to 3000 x g

STERILITY

- Gamma irradiation

Handling / Directions For Use

Use as directed in protocol of choice.

Notes and Tips

For chemical resistance or thermoplastic property information, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products available from STEMCELL Technologies, visit www.stemcell.com or contact us at techsupport@stemcell.com.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2019 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Falcon is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.