

# Round-Bottom Polypropylene Tubes with Caps, 5 mL

## Polypropylene tubes with caps

Catalog # 100-0089  
38057

25 Tubes  
500 Tubes



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

Falcon® Round-Bottom Polypropylene Tubes, 5 mL are sterile, high-quality tubes ideal for a wide variety of sample preparation, sample storage, and sample centrifugation tasks. They are required for select protocols in cell isolation, flow cytometry, and fluorescence-activated cell sorting (FACS). Tubes are RNase-/DNase-free.

Catalog #100-0089, 5 mL - 25 Tubes with caps; 25 Tubes/pack (Corning Catalog #352063)

Catalog #38057, 5 mL - 500 Tubes with caps; 25 Tubes/pack (Corning Catalog #352063)

### APPLICATIONS

- Cell separation
- Cell culture
- Sample storage
- Liquid handling

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** For Catalog #100-0089: Stable until expiry date (EXP) on STEMCELL label.

For Catalog #38057: Stable until expiry date on manufacturer label.

## Features

### MATERIALS

- Tube: clear polypropylene
- Cap: polyethylene
- Dual-position snap cap with a tight seal or a loose position for sterile aerobic culture
- Tube and cap are biologically inert, RNase-/DNase-free, and non-pyrogenic

### DIMENSIONS

- 12 x 75 mm (outer diameter x height)

### NOMINAL VOLUME

- 5 mL
- No graduations

### CENTRIFUGAL FORCE

- Up to 3000 x g

### STERILITY

- Gamma irradiation

## Handling / Directions For Use

Use as directed in protocol of choice.

## Notes and Tips

For chemical resistance or thermoplastic property information, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products available from STEMCELL Technologies, visit [www.stemcell.com](http://www.stemcell.com) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2019 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Falcon is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.