

Small Molecules

Hydrocortisone

Enhances proliferation and differentiation of endothelial, epithelial, mesenchymal, and oligodendrocyte cells

Catalog #74142
74144

100 mg
1 g



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Hydrocortisone is a glucocorticoid hormone with anti-inflammatory and immunosuppressive effects (Kang et al.; Langhoff & Ladefoged). It is also commonly used as a supplement for endothelial, epithelial, mesenchymal, or oligodendrocyte cell culture media, as it supports growth and differentiation. This product is supplied as the sodium salt of hydrocortisone 21-hemisuccinate, which is a more water-soluble form of hydrocortisone. Hydrocortisone is required as a supplement in MyeloCult™ H5100 (Human; Catalog #05100) for long-term cultures and the long-term culture-initiating cell (LTC-IC) assay. For complete instructions, refer to the Technical Manual for Human (Document #28412) or Mouse (Document #10000005500) LTC-IC Assays, available at www.stemcell.com, or contact us to request a copy.

Molecular Name:	Hydrocortisone
Alternative Names:	Cortisol 21-hemisuccinate; Hydrocortisone sodium succinate; ST51037292; U-4905
CAS Number:	125-04-2
Chemical Formula:	C ₂₅ H ₃₃ O ₈ • Na
Molecular Weight:	484.5 g/mol
Purity:	≥ 98%
Chemical Name:	21-(3-carboxy-1-oxopropoxy)-11β,17-dihydroxy-pregn-4-ene-3,20-dione, monosodium salt
Structure:	



Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at -20°C as supplied. Protect from prolonged exposure to light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
Solubility:	<ul style="list-style-type: none">• PBS (pH 7.2) ≤ 20 mM• DMSO ≤ 60 mM• Absolute ethanol ≤ 40 mM <p>For example, to prepare a 10 mM stock solution in DMSO, resuspend 10 mg in 2.06 mL of DMSO.</p> <p>Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.</p> <p>For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.</p>

Published Applications

MAINTENANCE AND SELF-RENEWAL

- Enhances proliferation and self-renewal of limbal stem cells in explant cultures (Yu et al.).
- Supports long-term growth of immature human myeloid cells in culture (Salahuddin et al.).
- Increases the activity of glucose 6-phosphate dehydrogenase in mouse mammary epithelial cells in culture in the presence of insulin (Oka & Perry).

DIFFERENTIATION

- With ATRA, EGF, KGF, and HGF, promotes epithelial differentiation of adipose-derived stem cells (Li et al.).
- Promotes differentiation of mouse embryonic stem (ES) cell-derived definitive endoderm toward lung alveolar epithelial cells (Mokhber Dezfouli et al.).

IMMUNOLOGY

- Inhibits the bioactivity of the proinflammatory cytokine interleukin 6 ($IC_{50} = 6.7 \text{ mM}$; Kang et al.).
- Suppresses the activity on phytohemagglutinin (PHA)-stimulated lymphocytes in vitro (Langhoff & Ladefoged).

References

- Kang BS et al. (2001) Inhibitory effects of anti-inflammatory drugs on interleukin-6 bioactivity. *Biol Pharm Bull* 24(6): 701–3.
- Langhoff E & Ladefoged J. (1983) Relative immunosuppressive potency of various corticosteroids measured in vitro. *Eur J Clin Pharmacol* 25(4): 459–62.
- Li H et al. (2012) Effects of multiple agents on epithelial differentiation of rabbit adipose-derived stem cells in 3D culture. *Tissue Eng Part A* 18(17–18): 1760–70.
- Mokhber Dezfouli MR et al. (2019) Hydrocortisone promotes differentiation of mouse embryonic stem cell-derived definitive endoderm toward lung alveolar epithelial cells. *Cell J* 20(4): 469–76.
- Oka T & Perry JW. (1974) Studies on the function of glucocorticoid in mouse mammary epithelial cell differentiation in vitro. Stimulation of glucose 6-phosphate dehydrogenase. *J Biol Chem* 249(11): 3586–91.
- Salahuddin SZ et al. (1981) Long-term suspension cultures of human cord blood myeloid cells. *Blood* 58(5): 931–8.
- Yu M et al. (2016) An important role for adenine, cholera toxin, hydrocortisone and triiodothyronine in the proliferation, self-renewal and differentiation of limbal stem cells in vitro. *Exp Eye Res* 152: 113–22.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules, or contact us at techsupport@stemcell.com.

This product is hazardous. Please refer to the Safety Data Sheet (SDS).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2022 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.