

# Cytokines

## Human Recombinant NT-4



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

### Neurotrophin 4

Catalog # 78093  
78093.1

10 µg  
100 µg

## Product Description

Neurotrophin-4 (NT-4) is a member of the nerve growth factor family which includes neurotrophin-3 (NT-3), brain-derived neurotrophic factor (BDNF), and nerve growth factor (NGF), all of which promote the differentiation, growth, and survival of peripheral and central nervous system neurons (Eide et al.). NT-4 binds and activates tropomyosin receptor kinase B (TrkB) at the cell surface; in doing so, it acts as a survival factor for certain populations of sensory neurons (Berkemeier et al.; Skaper). It has been shown that NT-4, together with BDNF, promotes neurite extension and maturation, as well as maintenance of differentiated cerebellar granule cells (Gao et al.).

## Product Information

**Alternative Names:** GLC10, GLC10, Neurotrophic factor 4, Neurotrophic factor 5, Neurotrophin-5, NT-4/5, NT-5, NTF4, NTF5  
**Accession Number:** P34130  
**Amino Acid Sequence:** MGVSETAPAS RRGELAVCDA VSGWVTD RRT AVDLRGREVE VLGEVPAAGG SPLRQYFFET RCKADNAEEG GPGAGGGGCR GVDRRHVWSE CKAKQSYVRA LTADAQGRVG WRWIRIDTAC VCTLLSRTGR A  
**Predicted Molecular Mass:** 14.0 kDa monomer; 28.1 kDa dimer  
**Species:** Human  
**Formulation:** Lyophilized from a sterile-filtered aqueous solution containing 0.1% trifluoroacetic acid.  
**Source:** *E. coli*

## Specifications


**Activity:** **Lot #1000039464 or higher:** The specific activity is  $\geq 5.0 \times 10^4$  units/mg ( $EC_{50} \leq 20$  ng/mL) as determined by a cell proliferation assay using a neuroblastoma cell line stably expressing TrkB (BR6).  
**All other lots:** The specific activity is  $\geq 3.3 \times 10^2$  units/mg ( $EC_{50} \leq 3.0$  µg/mL) as determined by a cell proliferation assay using C6 cells.  
**Purity:**  $\geq 95\%$   
**Endotoxin Level:** Measured by kinetic Limulus amoebocyte lysate (LAL) analysis and is  $\leq 1$  EU/µg protein.

## Preparation and Storage

**Storage:** Store at -20°C to -80°C.  
**Stability:** Stable as supplied for 12 months from date of receipt.  
**Preparation:** Centrifuge vial before opening. Reconstitute the product in sterile water to at least 0.1 mg/mL by pipetting the solution down the sides of the vial. Do not vortex.

**OPTIONAL:** After reconstitution, if product will not be used immediately, dilute with concentrated bovine serum albumin (BSA) to a final BSA concentration of 0.1%. The effect of storage of stock solution on product performance should be tested for each application. As a general guide, do not store at 2 - 8°C for more than 1 month or at -80°C for more than 3 months. Avoid repeated freeze-thaw cycles.

## Data



(A) The biological activity of Human Recombinant NT-4 was tested by its ability to promote the proliferation of a neuroblastoma cell line stably expressing TrkB (BR6). Cell proliferation was measured using a colorimetric assay method. The EC<sub>50</sub> is defined as the effective concentration of the growth factor at which cell proliferation is at 50% of maximum. The EC<sub>50</sub> in the above example is 5.6 ng/mL. For a representative bioactivity plot of C6 cell proliferation, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

(B) 1 µg of Human Recombinant NT-4 was resolved with SDS-PAGE under reducing (+) and non-reducing (-) conditions and visualized by Coomassie Blue staining. Human Recombinant NT-4 is a homodimer of 14 kDa subunits with a predicted total molecular mass of 28.1 kDa.

## Related Products

For a complete list of cytokines, as well as related products available from STEMCELL Technologies, visit [www.stemcell.com/cytokines](http://www.stemcell.com/cytokines) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## References

- Berkemeier LR et al. (1991) Neurotrophin-5: a novel neurotrophic factor that activates trk and trkB. *Neuron* 7(5): 857–66.
- Eide FF et al. (1993) Neurotrophins and their receptors--current concepts and implications for neurologic disease. *Exp Neurol* 121(2): 200–14.
- Gao WQ et al. (1995) Neurotrophin-4/5 (NT-4/5) and brain-derived neurotrophic factor (BDNF) act at later stages of cerebellar granule cell differentiation. *J Neurosci* 15(4): 2656–67.
- Skaper SD. (2008) The biology of neurotrophins, signalling pathways, and functional peptide mimetics of neurotrophins and their receptors. *CNS Neurol Disord Drug Targets* 7(1): 46–62.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2021 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.