

Human Recombinant IFN-beta (HEK293-expressed)

Interferon-beta

Catalog #78113	5 µg
Catalog #78113.1	25 µg

Product Description

Interferon beta (IFN-β), also known as type I interferon, is produced by fibroblasts, and in smaller amounts by plasmacytoid dendritic cells. Macrophages and endothelial cells secrete IFN-β in response to a viral infection (Reder & Feng). IFN-β binds to a receptor complex composed of IFNAR1 and IFNAR2, and initiates signal transduction via the JAK/STAT pathway; this culminates in the transcription and activation of many genes that control dendritic cell activation, T cell survival, NK cell activation, chemokine expression, lymph node retention, and antiproliferative and antiviral effects (Dunn et al.). IFN-β is a first-line treatment for multiple sclerosis. It suppresses Th17 cells by affecting expression of IL-4, IL-10, and IL-27. IFN-β was also shown to expand regulatory T cells and limit T cell trafficking to the central nervous system (Inoue & Shinohara). Out of the two variants of IFN-β (IFN-β1 and IFN-β3), this product is the IFN-β1 form.

Product Information

Alternative Names:	B cell interferon, Fibroblast interferon, IFNB1, Leukocyte interferon, Type I interferon
Accession Number:	P01574
Amino Acid Sequence:	MSYNLLGFLQ RSSNFQCQKL LWQLNGRLEY CLKDRMNFDI PEEIKQLQQF QKEDAALTIY EMLQNIFAIF RQDSSSTGWN ETIVENLLAN VYHQINHLKT VLEEKLEKED FTRGKLMSSL HLKRYYGRIL HYLKAKEYSH CAWTIVRVEI LRNFYFINRL TGYLRN
Predicted Molecular Mass:	20 kDa
Species:	Human
Product Formulation:	Lyophilized after dialysis against phosphate-buffered saline.
Source:	HEK293
Purity:	≥ 95%

Specifications

Activity:	The specific activity is $\geq 1.0 \times 10^7$ units/mg ($EC_{50} \leq 0.1$ ng/mL), as determined by a cell proliferation assay using TF-1 cells.
Endotoxin Level:	Measured by kinetic Limulus amoebocyte lysate (LAL) analysis and is ≤ 0.2 EU/µg protein.

Preparation and Storage

Stability and Storage:


Store at -80°C. Stable as supplied for 12 months from date of receipt.

Preparation:

Centrifuge vial before opening. Reconstitute the product in sterile water to at least 0.1 mg/mL by pipetting the solution down the sides of the vial. Do not vortex.

OPTIONAL: After reconstitution, if product will not be used immediately, dilute with concentrated bovine serum albumin (BSA) to a final BSA concentration of 0.1%. The effect of storage of stock solution on product performance should be tested for each application. As a general guide, do not store at 2 - 8°C for more than 1 week or at -20°C for more than 3 months. Avoid repeated freeze-thaw cycles.

Data



(A) The biological activity of Human Recombinant IFN-beta was tested by its ability to inhibit the proliferation of TF-1 cells. The EC50 is defined as the effective concentration of the growth factor at which cell proliferation is at 50% of maximum. The EC50 in the above example is less than 0.1 ng/mL. (B) 2 µg of Human Recombinant IFN-beta was resolved with SDS-PAGE under reducing (+) and non-reducing (-) conditions and visualized by Coomassie Blue staining. Human Recombinant IFN-beta has a predicted molecular mass of 20 kDa.

Related Products

For a complete list of cytokines or peptide pools, as well as related products available from STEMCELL Technologies, visit www.stemcell.com/cytokines or contact us at techsupport@stemcell.com.

References

Dunn GP et al. (2006) Interferons, immunity and cancer immunoediting. *Nat Rev Immunol* 6(11): 836–48.

Inoue M & Shinohara ML. (2013) The role of interferon-β in the treatment of multiple sclerosis and experimental autoimmune encephalomyelitis - in the perspective of inflammasomes. *Immunology* 139(1): 11–8.

Reder AT & Feng X. (2013) Aberrant type I interferon regulation in autoimmunity: opposite directions in MS and SLE, shaped by evolution and body ecology. *Front Immunol* 4: 281.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.