

# Cytokines

## Human Recombinant DKK-1

Dickkopf-related protein 1



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog #78208

#78208.1

#78208.2

10 µg

50 µg

1000 µg

## Product Description

Dickkopf-related protein 1 (DKK-1) is a member of the Dickkopf family and is a secreted protein that inhibits the canonical WNT pathway by competitive binding to low-density lipoprotein receptors (LRP)-5 and -6 with high affinity, thereby decreasing  $\beta$ -catenin protein stability (Niehrs). DKK-1 regulates embryonic development and contains two conserved cysteine-rich domains separated by a linker region and an N-terminal signal peptide (Krupnik et al.; Lieven et al.). It contains a His-residue tag at the carboxyl end of the polypeptide chain. A family of human DKK-related genes composed of DKK-1, DKK-2, DKK-3, and DKK-4 have been characterized together with a unique DKK-3 related protein termed Soggy (Krupnik et al.). DKK-1 has been shown to support the generation of myeloid-derived suppressor cells (MDSCs) and thus is a negative regulator of antitumor immune responses (D'Amico et al.). DKK-1 from thrombocytes is an important regulator of leukocyte infiltration and induces Th2 cell polarization and potentiates Th2 cell cytokine expression (Chae et al.). DKK-1 has also been shown to drive cardiac and retinal differentiation from induced pluripotent stem (iPS) cells (Lian et al.).

## Product Information

**Alternative Names:** Dickkopf-1, Dickkopf WNT signaling pathway inhibitor 1, hDkk-1, SK

**Accession Number:** O94907

**Amino Acid Sequence:** TLNSVLNSNA IKNLPPPLGG AAGHPGSAVS AAPGILYPPG NKYQTIDNYQ PYPCAEDEEC GTDEYCASPT RGGDAGVQIC LACRKRKRRC MRHAMCCPGN YCKNGICVSS DQNHFRGEIE ETITESFGND HSTLDGYSRR TTLSSKMYHT KGQEGSVCLR SSDCASGLCC ARHFWSKICK PVLKEGQVCT KHRRKGGSHGL EIFQRCYCGE GLSCRIQKDH HQASNSSLRH TCQRHHHHHH H

**Predicted Molecular Mass:** 26.6 kDa

**Species:** Human

**Formulation:** Lyophilized from a sterile-filtered solution containing phosphate-buffered saline.

**Source:** HEK293 cells

## Specifications

**Activity:** The specific activity is  $\geq 250$  units/mg ( $EC_{50} \leq 4$  µg/mL), as determined by the ability to inhibit Wnt3a-induced alkaline phosphatase production in C3H10T1/2 cells.

**Purity:**  $\geq 95\%$

**Endotoxin Level:** Measured by kinetic Limulus amoebocyte lysate (LAL) analysis and is  $\leq 0.2$  EU/µg protein.

## Preparation and Storage


**Storage:** Store at  $-80^{\circ}\text{C}$ .

**Stability:** Stable as supplied for 12 months from date of receipt.

**Preparation:** Centrifuge vial before opening. Reconstitute the product in sterile water to at least 0.1 mg/mL by pipetting the solution down the sides of the vial. Do not vortex.

OPTIONAL: After reconstitution, if product will not be used immediately, dilute with concentrated bovine serum albumin (BSA) to a final BSA concentration of 0.1%. The effect of storage of stock solution on product performance should be tested for each application. As a general guide, do not store at  $2 - 8^{\circ}\text{C}$  for more than 1 week or at  $-20^{\circ}\text{C}$  for more than 3 months. Avoid repeated freeze-thaw cycles

## Data



(A) The biological activity of Human Recombinant DKK-1 was tested by its ability to inhibit Wnt3a induced alkaline phosphatase production by C3H10T1/2 cells. Alkaline phosphatase production was measured using a fluorometric assay method. The EC<sub>50</sub> is defined as the effective concentration of the growth factor at which alkaline phosphatase production is at 50% of maximum. The EC<sub>50</sub> in the above example is less than 4  $\mu\text{g/mL}$ .

(B) 2  $\mu\text{g}$  of Human Recombinant DKK-1 was resolved with SDS-PAGE under reducing (+) and non-reducing (-) conditions and visualized by Coomassie Blue staining. Human Recombinant DKK-1 has a predicted molecular mass of 26.6 kDa and an apparent molecular mass of 40 kDa.

## Related Products

For a complete list of cytokines, as well as related products available from STEMCELL Technologies, visit [www.stemcell.com/cytokines](http://www.stemcell.com/cytokines), or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## References

- Chae W-J et al. (2016) The Wnt antagonist Dickkopf-1 promotes pathological type 2 cell-mediated inflammation. *Immunity* 44(2): 246–58.
- D'Amico L et al. (2016) Dickkopf-related protein 1 (Dkk1) regulates the accumulation and function of myeloid derived suppressor cells in cancer. *J Exp Med* 213(5): 827–40.
- Krupnik VE et al. (1999) Functional and structural diversity of the human Dickkopf gene family. *Gene* 238(2): 301–13.
- Lian X et al. (2012) Robust cardiomyocyte differentiation from human pluripotent stem cells via temporal modulation of canonical Wnt signaling. *Proc Natl Acad Sci USA* 109(27): E1848–57.
- Lieven O et al. (2010) The regulation of Dkk1 expression during embryonic development. *Dev Biol* 340(2): 256–68.
- Niehrs C. (2006) Function and biological roles of the Dickkopf family of Wnt modulators. *Oncogene* 25(57): 7469–81.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.