

Double Chamber Slide Kit



For the culture, fixation, and staining of collagen-based cultures

Catalog #04963 1 Kit
#04813 16/Box
#04911 48/Pack

Scientists Helping Scientists™ | www.stemcell.com

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

The Double Chamber Slide Kit (Catalog #04963) includes MegaCult™-C Double Chamber Slides (Catalog #04813) and Filter Cards and Spacers (Catalog #04911). Double chamber slides are recommended for performing colony assays in collagen gels using MegaCult™-C or other collagen-based media for the detection of hematopoietic progenitor cells. The chamber slide system allows for culture, fixation, and immunocytochemical staining of hematopoietic colony assays to be performed on the same slide. Filter cards and spacers are used for the dehydration of collagen-based gels following culture. The kit provides sufficient slides, filter cards, and spacers to perform 24 duplicate assays (two slides) or 48 single assays (one slide).

Product Information

| PRODUCT NAME | CATALOG # | SIZE | STORAGE | SHELF LIFE |
|--|-----------|---------|---------------------|----------------|
| Double Chamber Slide Kit <ul style="list-style-type: none">MegaCult™-C Double Chamber Slides (3 x 16 slides)Filter Cards and Spacers (48 of each) | 04963 | 1 Kit | Store at 15 - 25°C. | Not applicable |
| MegaCult™-C Double Chamber Slides | 04813 | 16/Box | Store at 15 - 25°C. | Not applicable |
| Filter Cards and Spacers | 04911 | 48/Pack | Store at 15 - 25°C. | Not applicable |

Directions for Use

For complete instructions for collagen-based colony assays of human megakaryocyte progenitor cells, refer to the Technical Manual: MegaCult™-C Assays for Quantitation of Human & Mouse Megakaryocytic Progenitor Cells (Document #10000005592), available at www.stemcell.com, or contact us to request a copy.

References

Dobo I et al. (1995) Collagen matrix: an attractive alternative to agar and methylcellulose for the culture of hematopoietic progenitors in autologous transplantation products. *J Hematother* 4(4): 281–7.

Mossuz P et al. (1998) Use of collagen for standardization of PBSC graft quality evaluation: a multicenter comparative analysis of commercial collagen-based and methylcellulose-based colony-forming unit (CFU) assay kits. *J Hematother* 7(4): 351–9.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and MegaCult are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.