

## Small Molecules

FK-866

Nicotinamide phosphoribosyltransferase inhibitor

Catalog #100-0263  
100-0264

5 mg  
10 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

FK-866 is a highly specific noncompetitive inhibitor of nicotinamide phosphoribosyltransferase (NAMPT), an enzyme that regulates NAD<sup>+</sup> biosynthesis from the natural precursor nicotinamide (Cameron et al.; Hasmann & Schemainda). In hepatocarcinoma cells, FK-866 activates AMP-activated protein kinase (AMPK) and downregulates mammalian target of rapamycin (mTOR) signaling (Schuster et al.).

Molecular Name:	FK-866
Alternative Names:	K 22.175
CAS Number:	658084-64-1
Chemical Formula:	C <sub>24</sub> H <sub>29</sub> N <sub>3</sub> O <sub>2</sub>
Molecular Weight:	391.5 g/mol
Purity:	≥ 98%
Chemical Name:	N-[4-(1-benzoyl-4-piperidinyl)butyl]-3-(3-pyridinyl)-2E-propenamide
Structure:	



## Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at -20°C as supplied. Protect from prolonged exposure to light. For long-term storage store with a desiccant. Stable as supplied for 12 months from date of receipt.
Solubility:	<ul style="list-style-type: none"><li>· PBS ≤ 1.2 mM</li><li>· DMSO ≤ 60 mM</li><li>· Absolute ethanol ≤ 100 mM</li></ul> For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 255 μL of DMSO.

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

## Published Applications

### METABOLISM

· FK-866 can be used to deplete NAD<sup>+</sup>, a central metabolism cofactor (Cameron et al.; Hasmann & Schemainda; Jadeja et al.).

### CANCER RESEARCH

· Depletes NAD<sup>+</sup> and induces delayed cell death by apoptosis in HepG2 human liver carcinoma cells (IC<sub>50</sub> = ~1 nM; Hasmann & Schemainda).

· Triggers dose-dependent cytotoxicity in multiple myeloma cells (Cea et al.).

· Induces autophagic death in neuroblastoma SH-SY5Y cells (IC<sub>50</sub> = 0.93 nM; Billington et al.).

## References

Billington RA et al. (2008) NAD depletion by FK866 induces autophagy. *Autophagy* 4(3): 385–7.

Cameron AM et al. (2019) Inflammatory macrophage dependence on NAD<sup>+</sup> salvage is a consequence of reactive oxygen species-mediated DNA damage. *Nat Immunol* 20: 420–32.

Cea M et al. (2012) Targeting NAD<sup>+</sup> salvage pathway induces autophagy in multiple myeloma cells via mTORC1 and extracellular signal-regulated kinase (ERK1/2) inhibition. *Blood* 120(17): 3519–29.

Hasmann M & Schemainda I. (2003) FK866, a highly specific noncompetitive inhibitor of nicotinamide phosphoribosyltransferase, represents a novel mechanism for induction of tumor cell apoptosis. *Cancer Res* 63(21): 7436–42.

Jadeja RN et al. (2018) Loss of NAMPT in aging retinal pigment epithelium reduces NAD<sup>+</sup> availability and promotes cellular senescence. *Aging* 10(6): 1306–23.

Schuster S et al. (2015) FK866-induced NAMPT inhibition activates AMPK and downregulates mTOR signaling in hepatocarcinoma cells. *Biochem Biophys Res Commun* 458(2): 334–40.

## Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit [www.stemcell.com/smallmolecules](http://www.stemcell.com/smallmolecules) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

**This product is hazardous. Please refer to the Safety Data Sheet (SDS).**

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2020 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.