

# Bovine Serum Albumin

**Bovine serum albumin for cell culture and biochemical applications**



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

STEMCELL bovine serum albumin (BSA) is highly purified serum albumin protein derived from cows and isolated from the fifth fraction of the Edwin Cohn purification methodology. BSA has numerous biochemical applications and is often used as a component of culture media. STEMCELL BSA is bovine virus free, protease free, and has low endotoxin levels.

## Ordering Information

CATALOG #	ORIGIN	SIZE
100-0174	AUS	10 kg
100-0175	AUS	1 kg
100-0181	AUS	100 g
100-0176	NZ	10 kg
100-0177	NZ	1 kg
100-0178	NZ	100 g

## Properties

**Storage:** Store at 2 - 8°C.  
**Shelf Life:** Stable as supplied until expiry date (EXP) on label.  
**Physical Appearance:** White to off-white amber lyophilized powder  
**CAS Number:** 9048-46-8  
**Predicted Molecular Mass:** ~67 kDa

## Directions For Use

Use as directed in the protocol of choice.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright ©2020 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.