

# CellAdhere™ Collagen I-Coated, 6-Well Flat-Bottom Plate



Scientists Helping Scientists™ | [www.stemcell.com](http://www.stemcell.com)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:info@stemcell.com) • [TECHSUPPORT@STEMCELL.COM](mailto:techsupport@stemcell.com)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

**Sterile flat-bottom, tissue culture-treated multiwell plate with lid**

Catalog # 100-0362	5 Plates
Catalog # 100-1054	10 x 5 Plates

## Product Description

CellAdhere™ Collagen I-Coated, 6-Well Flat-Bottom Plate is a sterile, clear polystyrene plate coated with high-purity collagen I derived from rat tail tendon for optimal cell attachment. These plates can be used with a variety of cell types and cell culture applications.

Catalog #100-0362, CellAdhere™ Collagen I-Coated, 6-Well Flat-Bottom Plate; 5 Plates/pack

Catalog #100-1054, CellAdhere™ Collagen I-Coated, 6-Well Flat-Bottom Plate; 10 x 5 Plates/pack

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** Stable until expiry (EXP) on label.

## Features

### MATERIALS

- Plate and lid: polystyrene
- Non-pyrogenic

### DIMENSIONS

- 127.8 x 85.2 x 22.4 mm (length x width x depth)

### CAPACITY

- Recommended working volume of 4 mL
- Approximate cell growth area of 53.2 cm<sup>2</sup>

### STERILITY

- Gamma irradiation

## Directions For Use

Use as directed in protocol of choice.

## Notes and Tips

For chemical resistance or thermoplastic property information, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products available from STEMCELL Technologies, visit [www.stemcell.com](http://www.stemcell.com) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2022 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and CellAdhere are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.