

CellAdhere™ Collagen I-Coated, 96-Well Flat-Bottom Plate



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Sterile flat-bottom, tissue culture-treated multiwell plate with lid

Catalog # 100-0366

5 Plates

Product Description

CellAdhere™ Collagen I-Coated, Flat-Bottom Plate is a sterile, clear polystyrene plate coated with high-purity collagen I derived from rat tail tendon for optimal cell attachment.

Catalog #100-0362, CellAdhere™ Collagen I-Coated, 6-Well Flat-Bottom Plate; 5 Plates/pack

Catalog #100-0363, CellAdhere™ Collagen I-Coated, 12-Well Flat-Bottom Plate; 5 Plates/pack

Catalog #100-0364, CellAdhere™ Collagen I-Coated, 24-Well Flat-Bottom Plate; 5 Plates/pack

Catalog #100-0365, CellAdhere™ Collagen I-Coated, 48-Well Flat-Bottom Plate; 5 Plates/pack

Catalog #100-0366, CellAdhere™ Collagen I-Coated, 96-Well Flat-Bottom Plate; 5 Plates/pack

APPLICATIONS

- Cell culture

Properties

Storage: Store at 15 - 25°C.

Shelf Life: Stable until expiry (EXP) on label.

Features

MATERIALS

- Plate and lid: polystyrene
- Non-pyrogenic

DIMENSIONS

- 127.8 x 85.2 x 16.55 mm (length x width x depth)

CAPACITY

- Recommended working volume of 0.2 mL
- Approximate cell growth area of 27.84 cm²

STERILITY

- Gamma irradiation

Handling / Directions For Use

Use as directed in protocol of choice.

Notes and Tips

For chemical resistance or thermoplastic property information, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products available from STEMCELL Technologies, visit www.stemcell.com or contact us at techsupport@stemcell.com.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2020 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and CellAdhere are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.