

RoboSep™-C PBS/EDTA Buffer

Cell separation buffer

Catalog #100-0188

3 x 3 L



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

RoboSep™-C PBS/EDTA Buffer is recommended for use with RoboSep™-C (Catalog #100-0185) when performing cell separation protocols.

Product Information

The following components are sold as a complete kit (Catalog #100-0188) and are not available for individual sale.

COMPONENT NAME	COMPONENT #	SIZE	STORAGE	SHELF LIFE	DESCRIPTION
RoboSep™-C PBS/EDTA Buffer	300-0129	3 L	Store at 15 - 25°C.	Stable for 12 months from date of manufacture (MFG) on label.	Dulbecco's phosphate-buffered saline (D-PBS) with 1 mM EDTA.

Directions For Use

RoboSep™-C PBS/EDTA Buffer has been optimized for use with RoboSep™-C cell separation reagents. When using RoboSep™-C PBS/EDTA Buffer as recommended medium in cell separation protocols, it may require the addition of a suitable protein source such as fetal bovine serum (FBS). For more information, refer to the Product Information Sheet (PIS) for the applicable cell separation kit or contact us at techsupport@stemcell.com. For the installation instructions, refer to the RoboSep™-C User Reference Manual (Document #1000008334), available at www.stemcell.com or contact us to request a copy.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO WWW.STEMCELL.COM/COMPLIANCE.

Copyright © 2020 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, EasySep, and RoboSep are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.