

Cl-Amidine

Inhibits protein-arginine deiminase (PAD)

Catalog #100-0518 1 mg

Catalog #100-0519 5 mg

Product Description

Cl-Amidine is a protein-arginine deiminase (PAD) inhibitor that irreversibly inactivates four subtypes of PAD ($k_{\text{inact}}/K_i = 37,000/\text{PAD1}$, $1,200/\text{PAD2}$, $2,000/\text{PAD3}$, and $13,000/\text{PAD4}$ $\text{M}^{-1}\text{min}^{-1}$; Knuckley et al.; Luo et al.; Slack et al.) by modifying the enzyme's active site. This product is supplied as the hydrochloride salt of the molecule.

Molecular Name: Cl-Amidine (Hydrochloride)

Alternative Names: Not applicable

CAS Number (Model): 1373232-26-8


Chemical Formula: $\text{C}_{14}\text{H}_{19}\text{ClN}_4\text{O}_2 \cdot \text{HCl}$

Molecular Weight: 347.2 g/mol

Purity: $\geq 95\%$

Chemical Name: N-[(1S)-1-(aminocarbonyl)-4-[(2-chloro-1-iminoethyl)amino]butyl]-benzamide, monohydrochloride

Structure:



Properties

Product Format:	A crystalline solid
Stability and Storage:	Product stable at -20°C as supplied. As a precaution, STEMCELL recommends storing all small molecules away from direct light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
Preparation:	<p>Solubility:</p> <ul style="list-style-type: none">• PBS (pH 7.2) ≤ 8.6 mM• DMSO ≤ 140 mM• Absolute ethanol ≤ 55 mM <p>For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 288 µL of DMSO.</p> <p>Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.</p> <p>For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO or absolute ethanol concentration above 0.1% due to potential cell toxicity.</p>

Published Applications

IMMUNOLOGY

- Inhibits neutrophil extracellular traps formation in neutrophils (Knight et al.).
- Prevents hypercitrullination of histone H3 in neutrophils (Knight et al.).

CANCER RESEARCH

- Exhibits cytotoxic effects toward human leukemia (IC₅₀ = 0.25 µM), breast (IC₅₀ = 0.05 µM), and colon (IC₅₀ = 1 µM) cancer cell lines (Slack et al.).

References

Knight JS et al. (2013) Peptidylarginine deiminase inhibition is immunomodulatory and vasculoprotective in murine lupus. *J Clin Invest* 123(7): 2981–93.

Knuckley B et al. (2010) Substrate specificity and kinetic studies of PADs 1, 3, and 4 identify potent and selective inhibitors of protein arginine deiminase 3. *Biochemistry* 49(23): 4852–63.

Luo Y et al. (2006) Inhibitors and inactivators of protein arginine deiminase 4: functional and structural characterization. *Biochemistry* 45(39): 11727–36.

Slack JL et al. (2011) Protein arginine deiminase 4: a target for an epigenetic cancer therapy. *Cell Mol Life Sci* 68(4): 709–20.

Related Products

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

Warning

This product is hazardous. Please refer to the Safety Data Sheet (SDS).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.