

Small Molecules

Camostat

Inhibits transmembrane serine protease (TMPRSS2)

Catalog #100-0552
100-0553

10 mg
50 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM


FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Camostat is a serine protease inhibitor with antiviral properties (Coote et al.; Hoffmann et al.). It inhibits the transmembrane serine protease TMPRSS2, which is involved in viral infections such as COVID-19 (Hoffmann et al.). This product is supplied as the mesylate salt of the molecule.

Molecular Name:	Camostat (Mesylate)
Alternative Names:	Foipan; FOY 305
CAS Number:	59721-29-8
Chemical Formula:	$C_{20}H_{22}N_4O_5 \cdot CH_3SO_3H$
Molecular Weight:	494.5 g/mol
Purity:	≥ 98%
Chemical Name:	4-[[4-[(aminoiminomethyl)amino]benzoyl]oxy]-benzeneacetic acid, 2-(dimethylamino)-2-oxoethyl ester, monomethanesulfonate

Structure:



Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at -20°C as supplied. Protect product from prolonged exposure to light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
Solubility:	• DMSO ≤ 50 mM For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 202 µL of DMSO. Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application. Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

Published Applications

IMMUNOLOGY

· Inhibits MCP-1 and TNF-alpha production in monocytes (Gibo et al.).

DISEASE MODELING

· Inhibits the proliferation and activation of pancreatic stellate cells in rats (Emori et al.).

References

Coote K et al. (2009) Camostat attenuates airway epithelial sodium channel function in vivo through the inhibition of a channel-activating protease. *J Pharmacol Exp Ther* 329(2): 764–74.

Emori Y et al. (2005) Camostat, an oral trypsin inhibitor, reduces pancreatic fibrosis induced by repeated administration of a superoxide dismutase inhibitor in rats. *J Gastroenterol Hepatol* 20(6): 895–9.

Gibo J et al. (2005) Camostat mesilate attenuates pancreatic fibrosis via inhibition of monocytes and pancreatic stellate cells activity. *Lab Invest* 85(1): 75–89.

Hoffmann M et al. (2020) SARS-CoV-2 cell entry depends on ACE2 and TMPRSS2 and is blocked by a clinically proven protease inhibitor. *Cell* 181(2): 271–80.e8.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

This product is hazardous. Please refer to the Safety Data Sheet (SDS).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2021 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.