

## Anti-SARS-CoV-2 Spike Protein S1 Receptor-Binding Domain Antibody, Clone Covi-2 (Blocking/Recombinant)



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Antibodies

Human monoclonal IgG1 antibody against SARS-CoV-2 (2019-nCoV) S protein (HEK293-expressed recombinant)

Catalog #100-0584

100 µg      1 mg/mL

## Product Description

The Covi-2 antibody reacts with Spike protein S1 receptor-binding domain (RBD) expressed by SARS-associated coronavirus 2 (SARS-CoV-2/2019-nCoV). The Spike (S) protein is a type I transmembrane glycoprotein present on the surface of coronaviruses. Entry of CoV into host cells is mediated by the S protein, where it interacts with the cell-surface receptor angiotensin-converting enzyme 2 (ACE2). In humans, ACE2 is expressed in several organs and tissues, including intestinal and respiratory epithelium. The S protein has two subunits, S1 and S2, where S1 primarily consists of the 193-amino-acid RBD and the N-terminal domain (NTD). The S2 domain is responsible for membrane fusion. During CoV infection, the S protein is cleaved into the N-terminal S1 subunit and C-terminal S2 subunit by host proteases, transforming its conformation from the pre-fusion to the post-fusion state. The S protein has been shown to play a key role in the induction of neutralizing antibody and T cell responses, which may lead to protective immunity. The RBD binds to ACE2, while the function of the NTD is not well understood. The Covi-2 antibody blocks the interaction of the S1 subunit with the ACE2 receptor and binds to a different epitope than the Covi-1 antibody.

<b>Target Antigen Name:</b>	SARS-CoV-2 Spike Protein S1 Receptor-Binding Domain (RBD)
<b>Alternative Names:</b>	SARS-CoV-2 Spike Protein, SARS-CoV-2 S Protein, SARS-CoV-2 S1 Protein
<b>Gene ID:</b>	N/A (Uniprot: P0DTC2)
<b>Species Reactivity:</b>	SARS-CoV-2 (2019-nCoV) Spike Protein S1. <ul style="list-style-type: none"><li>• Blocks interaction of SARS-CoV-2 Spike protein S1 RBD with human ACE2</li><li>• Detects a different epitope than Anti-SARS-CoV-2 Spike protein S1 RBD antibody, clone Covi-1</li></ul>
<b>Host Species:</b>	Human (recombinant antibody expressed in HEK293 cells)
<b>Clonality:</b>	Monoclonal
<b>Clone:</b>	Covi-2
<b>Isotype:</b>	IgG1
<b>Immunogen:</b>	RBD of SARS-CoV-2 Spike protein S1 (amino acids 330 - 530)
<b>Conjugate:</b>	Unconjugated

## Applications


<b>Verified:</b>	ELISA (Blocking)
<b>Reported:</b>	WB (Native)

Abbreviations: CellSep: Cell separation; ChIP: Chromatin immunoprecipitation; FA: Functional assay; FACS: Fluorescence-activated cell sorting; FC: Flow cytometry; ICC: Immunocytochemistry; IF: Immunofluorescence microscopy; IHC: Immunohistochemistry; IP: Immunoprecipitation; RIA: Radioimmunoassay; WB: Western blotting

## Properties

<b>Formulation:</b>	Phosphate-buffered solution
<b>Purification:</b>	The antibody was purified by affinity chromatography.
<b>Stability and Storage:</b>	Product stable at -20°C when stored undiluted. Stable until expiry date (EXP) on label.
<b>Directions for Use:</b>	It is recommended that the antibody be titrated for optimal performance for each application.

### Data



(A) ELISA binding assay showing various concentrations (x-axis) of Anti-SARS-CoV-2 Spike Protein S1 Receptor Binding Domain Antibody, Clone Covi-2 (Blocking/Recombinant) (black line) or a negative control antibody (grey line) blocking biotinylated recombinant human ACE2 from binding to an ELISA plate coated with recombinant SARS-CoV-2 (2019-nCoV) Spike protein S1 RBD. The binding percentage of biotinylated recombinant human ACE2 to the Spike protein S1 RBD (y-axis) was detected using streptavidin horseradish peroxidase (HRP). (B) ELISA binding assay showing Anti-SARS-CoV-2 Spike Protein S1 Receptor-Binding Domain Antibody, Clone Covi-2 (Blocking/Recombinant) (black line) or a negative control antibody (grey line) binding to a recombinant SARS-CoV-2 (2019-nCoV) Spike protein S1 RBD.

### Related Products

For a complete list of antibodies, including other conjugates, sizes, and clones, as well as related products available from STEMCELL Technologies, visit [www.stemcell.com/antibodies](http://www.stemcell.com/antibodies) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2020 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.