

Small Molecules

Lactacystin

Inhibits proteasomes

Catalog #100-0885

0.1 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Lactacystin is a metabolite isolated from *Streptomyces* and is an irreversible proteasome inhibitor that binds to the catalytic β -subunit ($IC_{50} = 4.8 \mu\text{M}$; Ōmura & Crump). It also induces neurite outgrowth in mouse neuroblastoma and apoptosis in human monoblast U937 cells (Imajoh-Ohmi et al.; Ōmura & Crump). Lactacystin can cause Parkinson's-like symptoms when injected into the substantia nigra of mice, making it useful to model protein degradation defects in Parkinson's disease (Savolainen et al.).

Alternative Names:	Not applicable
CAS Number:	133343-34-7
Chemical Formula:	$C_{15}H_{24}N_2O_7S$
Molecular Weight:	376.4 g/mol
Purity:	$\geq 98\%$
Chemical Name:	3S-hydroxy-2R-(1-hydroxy-2-methylpropyl)-4R-methyl-5-oxo-2-pyrrolidinecarboxylate-N-acetyl-L-cysteine
Structure:	



Properties

Physical Appearance:	A clear film (not easily visible)
Storage:	Product stable at -20°C as supplied. Protect product from prolonged exposure to light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
Solubility:	<ul style="list-style-type: none">· PBS (pH 7.2) $\leq 5.3 \text{ mM}$· DMSO $\leq 50 \text{ mM}$· Absolute ethanol $\leq 2.6 \text{ mM}$ <p>For example, to prepare a 10 mM stock solution in DMSO, resuspend 0.1 mg in 27 μL of DMSO.</p> <p>Prepare stock solution fresh before use; add diluent directly to the vial and vortex vigorously. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.</p> <p>For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.</p>

Published Applications

DIFFERENTIATION

· Induces neuroblastoma cell differentiation (Ōmura et al.).

References

Imajoh-Ohmi S et al. (1995) Lactacystin, a specific inhibitor of the proteasome, induces apoptosis in human monoblast U937 cells. *Biochem Biophys Res Commun* 217(3): 1070–7.

Ōmura S et al. (1991) Structure of lactacystin, a new microbial metabolite which induces differentiation of neuroblastoma cells. *J Antibiot (Tokyo)* 44(1): 117–8.

Ōmura S & Crump A. (2019) Lactacystin: first-in-class proteasome inhibitor still excelling and an exemplar for future antibiotic research. *J Antibiot (Tokyo)* 72(4): 189–201.

Savolainen MH et al. (2017) Nigral injection of a proteasomal inhibitor, lactacystin, induces widespread glial cell activation and shows various phenotypes of Parkinson's disease in young and adult mouse. *Exp Brain Res* 235(7): 2189–202.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2021 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.