

Cytokines

Human Recombinant R-Spondin-3

R-Spondin-3, His tag



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog #100-0934

50 µg

Product Description

R-Spondin-3 is a member of thrombospondin type 1 repeat (TSR-1) superfamily that is involved in the canonical Wnt/ β -catenin signaling pathway (de Lau et al.). R-spondin proteins are characterized by two furin-like repeats at the amino terminus and thrombospondin domain located near the carboxyl terminus (de Lau et al.). R-spondin-3 expression is associated with ovarian cancer (Gu et al.), prostate cancer (Mesci et al.), and differentiation of intestinal epithelial cells in diabetes mellitus (Shan et al.). In a transgenic mice model, the expression of R-Spondin-3 induces the expansion of Lgr5+ stem cells, Paneth cells, and Lgr4+ cells, promoting the intestinal stem cell compartment (Hilkens et al.). This protein contains a His-residue tag at the carboxyl end of the polypeptide chain.

Product Information

Alternative Names: Cristin 1, hRspo3, PWTSR, R-Spondin 3, Roof plate-specific spondin-3, RSPO3, RSpondin 3, Thrombospondin type-1 domain-containing protein 2, THSD2

Accession Number: Q9BXY4-1 (Gln22-Ala147) was expressed with a polyhistidine tag at the C-terminus

Amino Acid Sequence: QNASRGRRQR RMHPNVSQGC QGGCATCSDY NGCLSCPKRL FFALERIGMK QIGVCLSSCP SGYYGTRYPD INKCTKCKAD CDTCFNKNFC TKCKSGFY LH LGKCLDNCE GLEANNHTME CVSIVAHHHH HHHHHH

Predicted Molecular Mass: 15.3 kDa

Species: Human

Formulation: Lyophilized from sterile PBS, pH 7.4. Trehalose (5% - 8%), mannitol, and 0.01% TWEEN® 80 are normally added as protectants before lyophilization.

Source: HEK293

Specifications

Activity: Binding ability was measured in a functional ELISA. Immobilized Human Recombinant R-Spondin-3 at 2000 ng/mL can bind human RNF43 hFc with a linear range of 3 - 12 ng/mL.

Purity: \geq 93%

Endotoxin Level: Measured by kinetic Limulus amoebocyte lysate (LAL) analysis and is \leq 1.0 EU/ μ g protein.


Preparation and Storage

Storage: Store at -20°C to -80°C.

Stability: Stable as supplied for 12 months from date of receipt.

Preparation: Centrifuge vial before opening. Reconstitute the product in sterile water to at least 0.1 mg/mL by pipetting the solution down the sides of the vial. Do not vortex. The effect of storage of stock solution on product performance should be tested for each application. As a general guide, do not store at 2 - 8°C for more than 1 month or at -80°C for more than 3 months. Avoid repeated freeze-thaw cycles.

Data



(A) The binding activity of Human Recombinant R-Spondin-3 was tested by functional ELISA with immobilized Human Recombinant R-Spondin-3 at 2000 ng/mL. Immobilized Human Recombinant R-Spondin-3 can bind human RNF43 hFc with a linear range of 3 - 12 ng/mL.

(B) Human Recombinant R-Spondin-3 was resolved with SDS-PAGE under reducing (+) conditions and visualized by Coomassie Blue staining. Human Recombinant R-Spondin-3 has a predicted molecular mass of 15.3 kDa, but the apparent molecular mass is approximately 22 kDa to 27 kDa due to different glycosylation.

Related Products

For a complete list of cytokines, as well as related products available from STEMCELL Technologies, visit www.stemcell.com/cytokines or contact us at techsupport@stemcell.com.

References

- Gu H et al. (2020) RSPO3 is a marker candidate for predicting tumor aggressiveness in ovarian cancer. *Ann Transl Med* 8(21): 1351.
- Hilkens J et al. (2017) RSPO3 expands intestinal stem cell and niche compartments and drives tumorigenesis. *Gut* 66(6): 1095–105.
- de Lau WBM et al. (2012) The R-spondin protein family. *Genome Biol* 13(242).
- Mesci A et al. (2019) RSPO3 is a prognostic biomarker and mediator of invasiveness in prostate cancer. *J Transl Med* 17(1): 1–11.
- Shan TD et al. (2021) RSPO3 regulates the abnormal differentiation of small intestinal epithelial cells in diabetic state. *Stem Cell Res Ther* 12(1): 1–11.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2021 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. TWEEN is a registered trademark of Croda Americas LLC. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.