

Dyes and Stains

FITC-C6-DEVD-FMK

For detection of active caspase-3 in mammalian cells undergoing apoptosis

Catalog #100-0924

100 µg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Caspases are a family of proteases that play a central role in cellular apoptosis through the cleavage of select proteins, and result in the disassembly of the cell (Thornberry & Lazebnik). This process is switched on by pro-apoptotic signals from death receptors or stimuli such as cytotoxic reagents, which activate the initiator caspases and lead to activation of the effector caspases—caspase-3 and caspase-7 (Thornberry & Lazebnik). Caspase-3 has substrate selectivity for the peptide sequence Asp-Glu-Val-Asp (DEVD). FITC-C6-DEVD-FMK is a non-toxic, cell-permeable fluorescent substrate that irreversibly binds to activated caspase-3 in apoptotic cells. The FITC-labeled substrate can be detected with high sensitivity using flow cytometry. FITC-C6-DEVD-FMK is useful for screening caspase-3 inhibitors and quantifying apoptotic cells containing active caspase-3.

Molecular Weight:	994.99 g/mol
Excitation Wavelength:	491 nm
Emission Wavelength:	516 nm
Extinction Coefficient:	73,000 cm ⁻¹ M ⁻¹
Correction Factor:	0.254
Structure:	



Properties

Storage:	Store at -20°C.
Shelf Life:	Stable until expiry date (EXP) on box label. Protect product from prolonged exposure to light.
Format:	Solid

Directions for Use

Please read the entire protocol before proceeding.

Preparation of FITC-C6-DEVD-FMK Working Solution

1. To prepare a stock solution, dissolve FITC-C6-DEVD-FMK in dimethyl sulfoxide (DMSO) at 2 - 5 mM.
NOTE: If not used immediately, aliquot and store at -20°C. After thawing aliquots, use immediately; do not re-freeze.
2. To prepare a FITC-C6-DEVD-FMK working solution, dilute the stock solution (prepared in step 1) to 20 µM with Hanks' Balanced Salt Solution with 20 mM HEPES (HHBS). Use the working solution immediately; do not store.

Staining Cells

1. To induce apoptosis, treat cells with test compounds and incubate at 37°C and 5% CO₂ for an appropriate length of time.
Example: Jurkat cells treated with camptothecin require 4 - 6 hours of incubation for apoptosis to occur.
NOTE: The optimal cell density to induce apoptosis should be determined for different cell lines.
2. Prepare positive and negative controls as follows:
 - a. Positive control: camptothecin or staurosporine at 5X or 10X concentration
NOTE: The optimal concentration should be determined for different cell lines.
 - b. Negative control 1: untreated cells
 - c. Negative control 2: untreated cells + equal volume of FITC-C6-DEVD-FMK working solution
3. Add an equal volume of the FITC-C6-DEVD-FMK working solution to the cell solution and incubate at 37°C and 5% CO₂ for at least 1 hour.
NOTE: The optimal incubation time should be determined for different cell lines.
4. Centrifuge cells at 100 - 150 x *g* for 5 minutes. Remove supernatant and wash cells with HHBS at least once.
NOTE: Centrifugation speed may need to be optimized to ensure apoptotic cells are sufficiently pelleted.
5. Monitor the fluorescence intensity using a flow cytometer or a fluorescence microscope.

References

Thornberry NA & Lazebnik Y. (1998) Caspases: enemies within. *Science* 281(5381): 1312–6.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2022 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.