

# Cell Scraper, Scraper Blade, 25 cm

## Cell scraper with ABS handle and TPE blade

Catalog #200-0592

100 Scrapers



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

Cell Scraper, Scraper Blade, 25 cm has a sturdy acrylonitrile butadiene styrene (ABS) handle and a flexible thermoplastic elastomer (TPE) blade that prevents damaging cells or scratching surfaces. This product is ideal for harvesting cells and for use with Tissue Culture-Treated Dishes (Catalog #38046) and Cell Culture Flasks (Catalog #200-0500/200-0501). The scraper blade is positioned parallel to the handle and scrapes with a windshield-wiper motion.

Catalog #200-0592, 25 cm - 100 Scrapers, individually wrapped

### APPLICATIONS

- Cell culture

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** Stable until expiry date on manufacturer label. Avoid direct exposure to sunlight and high atmospheric humidity as these conditions may directly impact the product's packaging over time.

## Features

### MATERIALS

- Handle: acrylonitrile butadiene styrene (ABS), meets USP, Class VI requirements
- Blade: thermoplastic elastomer (TPE), meets USP, Class VI requirements
- RNase/DNase-free and non-pyrogenic
- All material resins are qualified using USP and/or ISO 10993 standards for cytotoxicity and have been shown to be non-toxic

### DIMENSIONS

- Handle: 25 cm
- Scraper blade: 20 mm

### STERILITY

- Gamma irradiation

## Directions for Use

Use as directed in protocol of choice.

## Notes and Tips

For chemical resistance or thermoplastic property information, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products, visit [www.stemcell.com](http://www.stemcell.com), or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2022 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness.