

# Falcon® 60 mm Dish, Tissue Culture-Treated



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Sterile tissue culture-treated dish

Catalog #200-0623

500 Dishes

## Product Description

Falcon® 60 mm Dish, Tissue Culture-Treated is ideal for attachment-dependent cell and tissue culture and distortion-free microscopic imaging. These polystyrene dishes are transparent and feature a unique easy-grip design for easy handling and storage.

Catalog #200-0623, 60 mm - 500 Dishes; 20 Dishes/pack (Corning Catalog #353002)

### APPLICATIONS

- Cell culture
- Microscopy

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** Stable until expiry date on manufacturer label.

## Features

### MATERIALS

- Dish and lid: clear and transparent polystyrene
- Tissue culture-treated for optimal cell attachment
- Non-pyrogenic

### DIMENSIONS

- 54.81 x 13.26 mm (diameter x height)

### CAPACITY

- Recommended working volume of 6 - 7 mL
- Approximate cell growth area of 21.29 cm<sup>2</sup>

### STERILITY

- Irradiation

## Directions for Use

Use as directed in protocol of choice.

## Notes and Tips

For chemical resistance or thermoplastic property information, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products, visit [www.stemcell.com](http://www.stemcell.com), or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Falcon is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness.