

Small Molecules

Fisetin

Sirtuin 1 activator

Catalog #100-1126

1 g



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Fisetin is a flavonoid found in plants and vegetables with senolytic activities (Yousefzadeh et al.). Fisetin upregulates the expression of sirtuin 1 (SIRT1), which is a member of the sirtuins superfamily histone deacetylases (HDAC) that helps inhibit many transcription factors, such as p53, nuclear factor kappa B (NF- κ B), and forkhead box O (FOXO) (Ayissi et al.). Fisetin inhibits DNA methyltransferase 1 (DNMT1)-mediated DNA methylation (IC_{50} = 3.5 μ M) through direct inhibition of the DNA methylase (Lee et al.).

Alternative Names:	CI-75620; NSC 407010; NSC 656275
CAS Number:	528-48-3
Chemical Formula:	C ₁₅ H ₁₀ O ₆
Molecular Weight:	286.2 g/mol
Purity:	≥ 90%
Chemical Name:	2-(3,4-dihydroxyphenyl)-3,7-dihydroxy-4H-1-benzopyran-4-one
Structure:	



Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at -20°C as supplied. As a precaution, STEMCELL recommends storing all small molecules away from direct light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
Solubility:	<ul style="list-style-type: none">• DMSO ≤ 100 mM• Absolute ethanol ≤ 17 mM <p>For example, to prepare a 10 mM stock solution in DMSO, resuspend 10 mg in 3.49 mL of DMSO.</p> <p>Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported; however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.</p> <p>Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO or absolute ethanol concentration above 0.1% due to potential cell toxicity.</p>

Published Applications

MAINTENANCE AND SELF-RENEWAL

- Enhances survival of PC12 cells from tunicamycin (Tm)-mediated cell death by mitigating apoptosis, autophagy, and reactive oxygen species (ROS) production (Yen et al.).

IMMUNOLOGY

- Downregulates phorbol-12-myristate 13-acetate plus calcium ionophore A23187 (PMACI)-stimulated gene expression and production of tumor necrosis factor- α (TNF- α), interleukin (IL)-1 β , IL-4, IL-6, and IL-8 in human mast cells (Park et al.).

CANCER RESEARCH

- Inhibits the PI3K/Akt and mTOR pathway in prostate, lung, myeloma, and melanoma cancer cells (Adhami et al.; Syed et al.).
- Known to possess antioxidant activity. Mechanistically, fisetin has been shown to attenuate H₂O₂-induced cell damage by scavenging ROS and activating GSH pathways in vitro (Kang et al.; Zhang et al.).

DISEASE MODELING

- Induces ERK activation which provides neuroprotective activities in Huntington's disease models (Maher et al.).
- Reduces senescence and age-related pathologies in yeast, simple eukaryotes, and mice (Gryniewicz and Demchuk; Wood et al.; Yousefzadeh et al.).

References

- Adhami VM et al. (2012) Dietary flavonoid fisetin: a novel dual inhibitor of PI3K/Akt and mTOR for prostate cancer management. *Biochem Pharmacol* 84(10): 1277–81.
- Ayissi VBO et al. (2014) Epigenetic effects of natural polyphenols: A focus on SIRT1-mediated mechanisms. *Mol Nutr Food Res* 58(1): 22–32.
- Gryniewicz G and Demchuk OM. (2019) New perspectives for fisetin. *Front Chem* 7: 697.
- Kang KA et al. (2014) Fisetin attenuates hydrogen peroxide-induced cell damage by scavenging reactive oxygen species and activating protective functions of cellular glutathione system. *Vitr Cell Dev Biol - Anim* 50(1): 66–74.
- Lee WJ et al. (2005) Mechanisms for the inhibition of DNA methyltransferases by tea catechins and bioflavonoids. *Mol Pharmacol* 68(4): 1018–30.
- Maher P et al. (2011) ERK activation by the polyphenols fisetin and resveratrol provides neuroprotection in multiple models of Huntington's disease. *Hum Mol Genet* 20(2): 261–70.
- Park H-H et al. (2007) Anti-inflammatory activity of fisetin in human mast cells (HMC-1). *Pharmacol Res* 55(1): 31–7.
- Syed D et al. (2013) Inhibition of Akt/mTOR Signaling by the dietary flavonoid fisetin. *Anticancer Agents Med Chem* 13(7): 995–1001.
- Wood JG et al. (2004) Sirtuin activators mimic caloric restriction and delay ageing in metazoans. *Nature* 430(7000): 686–9.
- Yen J-H et al. (2017) Fisetin protects PC12 cells from tunicamycin-mediated cell death via reactive oxygen species scavenging and modulation of Nrf2-driven gene expression, SIRT1 and MAPK signaling in PC12 cells. *Int J Mol Sci* 18(4): 852.
- Yousefzadeh MJ et al. (2018) Fisetin is a senotherapeutic that extends health and lifespan. *EBioMedicine* 36: 18–28.
- Zhang H et al. (2019) Nrf2–ARE signaling acts as master pathway for the cellular antioxidant activity of fisetin. *Molecules* 24(4): 708.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules, or contact us at techsupport@stemcell.com.

This product is hazardous. Please refer to the Safety Data Sheet (SDS).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.