

Small Molecules

MK-4827 (Tosylate)

PARP1 and PARP2 inhibitor

Catalog #100-1160

10 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

MK-4827 (Tosylate) is a potent, selective poly ADP ribose polymerase (PARP) inhibitor with affinity for both PARP1 and PARP2 ($IC_{50} = 3.8$ and 2.1 nM; Jones et al.). MK-4827 demonstrates antiproliferative effects on cancer cells with BRCA1 or BRCA2 mutations (Yuan et al.). PARPs are a family of enzymes that repair DNA breaks through the base excision repair (BER); MK-4827 inhibits this activity, leading to DNA damage and cell death (Yuan et al.).

Alternative Names:	Not applicable
CAS Number:	1038915-73-9
Chemical Formula:	$C_{19}H_{20}N_4O \cdot C_7H_8O_3S$
Molecular Weight:	492.6 g/mol
Purity:	$\geq 98\%$
Chemical Name:	2-[4-(3S)-3-piperidinylphenyl]-2H-indazole-7-carboxamide, 4-methylbenzenesulfonate
Structure:	



Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at $-20^{\circ}C$ as supplied. As a precaution, STEMCELL recommends storing all small molecules away from direct light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
Solubility:	<ul style="list-style-type: none">• $DMSO \leq 60$ mM• Absolute ethanol ≤ 2 mM <p>For example, to prepare a 10 mM stock solution in DMSO, resuspend 10 mg in 2.03 mL of DMSO.</p> <p>Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported; however, as a general guide we recommend storage in DMSO at $-20^{\circ}C$. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.</p> <p>Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO or absolute ethanol concentration above 0.1% due to potential cell toxicity.</p>

Published Applications

CANCER RESEARCH

- Antiproliferative effect against MDA-MB-436 and HeLa cells expressing BRCA1 mutation (Jones et al., 2009).
- Inhibits growth in human MDA-MB-436 tumor cells with BRCA 1 mutant (Jones et al., 2015).

References

Jones P et al. (2009) Discovery of 2-{4-[(3 S)-Piperidin-3-yl]phenyl}-2 H -indazole-7-carboxamide (MK-4827): a novel oral poly (ADP-ribose) polymerase (PARP) inhibitor efficacious in BRCA-1 and -2 mutant tumors. J Med Chem 52(22): 7170–85.

Jones P et al. (2015) Niraparib: a poly (ADP-ribose) polymerase (PARP) inhibitor for the treatment of tumors with defective homologous recombination. J Med Chem 58(8): 3302–14.

Yuan Y et al. (2011) Novel targeted therapeutics: inhibitors of MDM2, ALK and PARP. J Hematol Oncol 4(1): 16.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules, or contact us at techsupport@stemcell.com.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.