

# CellAdhere™ Gelatin-Coated, 96-Well Flat-Bottom Plate



Scientists Helping Scientists™ | [www.stemcell.com](http://www.stemcell.com)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog #100-1211  
#100-1227

5 Plates  
10 x 5 Plates

## Product Description

CellAdhere™ Gelatin-Coated, 96-Well Flat-Bottom Plate is a sterile, optically clear polystyrene plate coated with porcine-derived gelatin (200 µg/cm<sup>2</sup>) for optimal cell attachment.

Catalog #100-1211, CellAdhere™ Gelatin-Coated, 96-Well Flat-Bottom Plate; 5 Plates/pack  
Catalog #100-1227, CellAdhere™ Gelatin-Coated, 96-Well Flat-Bottom Plate; 10 x 5 Plates/pack

### APPLICATIONS

- Cell culture

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** Stable until expiry date (EXP) on label.

## Features

### MATERIALS

- Coating: porcine-derived gelatin (200 µg/cm<sup>2</sup>)
- Plate and lid: clear and transparent polystyrene
- Non-pyrogenic

### DIMENSIONS

- Conforms to standard plate footprint and dimensions

### CAPACITY

- Recommended working volume of 0.075 - 0.2 mL per well
- Approximate cell growth area of 27.8 cm<sup>2</sup> per well

### STERILITY

- Gamma irradiation

## Directions for Use

Use as directed in protocol of choice.

## Notes and Tips

For chemical resistance or thermoplastic property information, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products, visit [www.stemcell.com](http://www.stemcell.com), or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and CellAdhere are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness.