

CellAdhere™ Gelatin-Coated Flask, 75 cm²

Sterile polystyrene T-75 cm² flask with vent cap

Catalog #100-1212
#100-1228

5 Flasks
10 x 5 Flasks



Scientists Helping Scientists™ | www.stemcell.com

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713
INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM
FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

CellAdhere™ Gelatin-Coated Flask, 75 cm² with vent cap is a sterile, optically clear polystyrene flask coated with porcine-derived gelatin (200 µg/cm²) for optimal cell attachment. The canted neck of the flask offers good pipette and cell scraper access to the corners.

Catalog #100-1212, CellAdhere™ Gelatin-Coated Flask, 75 cm²; 5 Flasks/pack
Catalog #100-1228, CellAdhere™ Gelatin-Coated Flask, 75 cm²; 10 x 5 Flasks/pack

APPLICATIONS

- Cell culture

Properties

Storage: Store at 15 - 25°C.

Shelf Life: Stable until expiry date (EXP) on label.

Features

MATERIALS

- Coating: porcine-derived gelatin (200 µg/cm²)
- Flask: clear and transparent polystyrene
- Cap: polyethylene
- Non-pyrogenic

DIMENSIONS

- 154 x 81 x 37 mm (length x width x depth)

CAPACITY

- Recommended working volume of 15 - 22.5 mL
- Approximate cell growth area of 74.2 cm²

STERILITY

- Gamma irradiation

Directions for Use

Use as directed in protocol of choice.

Notes and Tips

For chemical resistance or thermoplastic property information, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products, visit www.stemcell.com, or contact us at techsupport@stemcell.com.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and CellAdhere are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness.