

# CellAdhere™ Collagen I-Coated Flask, 175 cm<sup>2</sup>



Scientists Helping Scientists™ | [www.stemcell.com](http://www.stemcell.com)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:info@stemcell.com) • [TECHSUPPORT@STEMCELL.COM](mailto:techsupport@stemcell.com)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog #100-1199  
#100-1225

5 Flasks  
10 x 5 Flasks

## Product Description

CellAdhere™ Collagen I-Coated Flask, 175 cm<sup>2</sup> with vent cap is a sterile, optically clear polystyrene flask coated with high-purity collagen I derived from rat tail tendon (10 µg/cm<sup>2</sup>) for optimal cell attachment. The canted neck of the flask offers good pipette and cell scraper access to the corners.

Catalog #100-1199, CellAdhere™ Collagen I-Coated Flask, 175 cm<sup>2</sup>; 5 Flasks/pack

Catalog #100-1225, CellAdhere™ Collagen I-Coated Flask, 175 cm<sup>2</sup>; 10 x 5 Flasks/pack

### APPLICATIONS

- Cell culture

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** Stable until expiry date (EXP) on label.

## Features

### MATERIALS

- Coating: collagen I derived from rat tail tendon (10 µg/cm<sup>2</sup>)
- Flask: clear and transparent polystyrene
- Cap: polyethylene
- Non-pyrogenic

### DIMENSIONS

- 229 x 115 x 50 mm (length x width x depth)

### CAPACITY

- Recommended working volume of 35 - 52.5 mL
- Approximate cell growth area of 175 cm<sup>2</sup>

### STERILITY

- Gamma irradiation

## Directions for Use

Use as directed in protocol of choice.

## Notes and Tips

For chemical resistance or thermoplastic property information, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products, visit [www.stemcell.com](http://www.stemcell.com), or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and CellAdhere are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness.