

# Cytokines

## Human Recombinant Complement Factor D



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog #100-1330

50 µg

## Product Description

Complement factor D, also known as adipsin, is a serine protease belonging to the S1 peptidase family. This protein is secreted by adipocytes into circulating blood and is also expressed by macrophages and monocytes (White et al.). It contains a His-residue tag at the carboxyl end of the polypeptide chain. Complement factor D is a component of the alternative pathway of the complement system and a part of the innate immune system, and it plays a role in the defense against infection (Barratt and Weitz). In the initiation phase of the complement pathway, complement factor D cleaves complement factor B (bound to component C3) to produce a complex known as C3 convertase. During the amplification phase, complement factor D cleaves complement factor B (bound to component C3b) to produce the C3bBb convertase and is involved in the propagation of complement activation. Apart from its immunological role, complement factor D is involved in other physiological processes, like the efficient clearing of damaged cell debris by phagocytes following acute liver injury (Cresci et al.). Deficiency in complement factor D is linked to increased susceptibility to pathogens like *Neisseria meningitidis* (Biesma et al.).

## Product Information

Alternative Names:	Adipsin, ADN, DF, Factor D protein, PFD
Accession Number:	P00746 (Ile26-Ala253) was expressed with a polyhistidine tag at the C-terminus.
Amino Acid Sequence:	ILGGREAEAH ARPYMASVQL NGAHLCGGVL VAEQWVLSAA HCLEDAADGK VQVLLGAHSL SQPEPSKRLY DVLRAVPHPD SQPDTIDHDL LLLQLSEKAT LGPAVRPLPW QRVDRDVAPG TLCDVAGWGI VNHAGRRPDS LQHVLLPVL D RACNRRRTHH DGAITERLMC AESNRRDSCK GDSGGPLVCG GVLEGVVTSG SRVCGNRKKP GIYTRVASYA AWIDSVLAAH HHHHHHHHHH
Predicted Molecular Mass:	25.8 kDa
Species:	Human
Formulation:	Lyophilized from sterile PBS, pH 7.4. Trehalose (5 - 8%), mannitol, and 0.01% TWEEN® 80 are normally added as protectants before lyophilization.
Source:	HEK293 cells


## Specifications

Activity:	Not available
Purity:	≥ 94%
Endotoxin Level:	Measured by kinetic Limulus amoebocyte lysate (LAL) analysis and is ≤ 1 EU/µg protein.

## Preparation and Storage

Storage:	Store at -20 to -80°C.
Stability:	Stable as supplied for 12 months from date of receipt.
Preparation:	Centrifuge vial before opening. Reconstitute the product in sterile water to at least 0.25 mg/mL by pipetting the solution down the sides of the vial. Do not vortex. As a general guide, do not store at 2 - 8°C for more than 1 month or at -80°C for more than 3 months. Avoid repeated freeze-thaw cycles.

## Data



Human Recombinant Complement Factor D was resolved with SDS-PAGE under reducing (+) conditions and visualized by Coomassie Blue staining. Human Recombinant Complement Factor D has a predicted molecular mass of 25.8 kDa and an apparent molecular mass of 33 kDa due to glycosylation.

## Related Products

For a complete list of cytokines, as well as related products available from STEMCELL Technologies, visit [www.stemcell.com/cytokines](http://www.stemcell.com/cytokines), or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## References

- Barratt J & Weitz I. (2021) Complement factor D as a strategic target for regulating the alternative complement pathway. *Front Immunol* 12: 712572.
- Biesma DH et al. (2001) A family with complement factor D deficiency. *J Clin Invest* 108(2): 233–40.
- Cresci GA et al. (2015) Alternative complement pathway component factor D contributes to efficient clearance of tissue debris following acute CCl<sub>4</sub>-induced injury. *Mol Immunol* 64(1): 9–17.
- White RT et al. (1992) Human adipsin is identical to complement factor D and is expressed at high levels in adipose tissue. *J Biol Chem* 267(13): 9210–3.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. TWEEN is a registered trademark of Croda Americas LLC. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.