

Antibodies

Anti-Human Synaptophysin (SYP) Antibody, Clone 249

Rabbit monoclonal IgG antibody
against human Synaptophysin, unconjugated



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog #100-1345

100 µL

Product Description

This rabbit monoclonal antibody (clone 249) reacts with human synaptophysin, a 38 kDa integral transmembrane protein possessing four transmembrane domains and a major small synaptic vesicle protein. It interacts with a variety of different motor proteins and synaptic vesicles, including dynamin I, myosin V, vesicular proton pump V-ATPase, and adaptor protein 1 (AP-1). Synaptophysin is extensively involved in vesicular functions such as protein sorting and priming, as well as synapse formation, exocytosis, and endocytosis. It is also a potent marker of axonal damage; synaptophysin was discovered colocalized with amyloid precursor protein (APP), a marker of axonal damage in patients with multiple sclerosis. Synaptophysin deletion in knock-out mice has been correlated with reduced spatial learning and impaired object novelty recognition.

Target Antigen Name:	Synaptophysin
Alternative Names:	MRX96, MRXSYP, XLID96
Gene ID:	6855
Species Reactivity:	Human
Host Species:	Rabbit
Clonality:	Monoclonal
Clone:	249
Isotype:	IgG
Immunogen:	Synthetic peptide of human SYP C-terminus
Conjugate:	Unconjugated

Applications


Verified:	ICC/IF
Special Applications:	This antibody clone has been verified for labeling synaptophysin-positive puncta in human pluripotent stem cell (hPSC)-derived forebrain-type neurons generated with STEMdiff™ Forebrain Neuron Differentiation Kit (Catalog #08600) and STEMdiff™ Forebrain Neuron Maturation Kit (Catalog #08605).

Abbreviations: CellSep: Cell separation; ChIP: Chromatin immunoprecipitation; FA: Functional assay; FACS: Fluorescence-activated cell sorting; FC: Flow cytometry; ICC: Immunocytochemistry; IF: Immunofluorescence microscopy; IHC: Immunohistochemistry; IHC-P: Immunohistochemistry (paraffin-embedded); IP: Immunoprecipitation; RIA: Radioimmunoassay; WB: Western blotting

Properties

Formulation:	Phosphate-buffered solution
Purification:	The antibody was purified by affinity chromatography.
Stability and Storage:	Product stable at -20°C when stored undiluted. Avoid repeated freeze-thaw cycles. Stable until expiry date (EXP) on label.
Directions for Use:	The suggested use of this antibody is: ICC/IF, 1 - 2 µg/mL. It is recommended that the antibody be titrated for optimal performance for each application. For antibody concentration, refer to the lot-specific Certificate of Analysis at www.stemcell.com/coa .

Data



Forebrain-type neurons were generated from hPSC-derived neural progenitor cells using STEMdiff™ Forebrain Neuron Differentiation Kit for 7 days and subsequently matured for the following 5 weeks using STEMdiff™ Forebrain Neuron Maturation Kit, then fixed and labeled with Anti-Human Synaptophysin Antibody, Clone 249, followed by Goat Anti-Rabbit IgG (H+L) Antibody, Polyclonal, iFluor™ 488 (green) (Catalog #100-1082). Nuclei were counterstained with DAPI (blue). Synaptophysin is concentrated as discrete puncta distributed along axonal processes that were labeled with Anti-Beta-Tubulin III Antibody, Clone TUJ1 (red) (Catalog #60052). Inset shows forebrain-type neurons labeled with a rabbit IgG isotype control antibody, followed by Goat Anti-Rabbit IgG (H+L) Antibody, Polyclonal, iFluor™ 488 (with DAPI staining).

Related Products

For a complete list of antibodies, including other conjugates, sizes, and clones, as well as related products available from STEMCELL Technologies, visit www.stemcell.com/antibodies, or contact us at techsupport@stemcell.com.

References

1. Disdier C et al. (2017) Brain inflammation, blood brain barrier dysfunction and neuronal synaptophysin decrease after inhalation exposure to titanium dioxide nano-aerosol in aging rats. *Scientific Reports* 7(1) : 1–13.
2. Seo H G et al. (2010) Early motor balance and coordination training increased synaptophysin in subcortical regions of the ischemic rat brain. *Journal of Korean Medical Science* 25(11) : 1638–45.
3. Gudi Viktoria et al. (2017) Synaptophysin is a reliable marker for axonal damage. *Journal of Neuropathology & Experimental Neurology* 76(2): 109–25.
4. Sarnat H B et al. (2010) Synaptophysin immunoreactivity in the human hippocampus and neocortex from 6 to 41 weeks of gestation. *Journal of Neuropathology & Experimental Neurology* 69(3): 234–45.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and STEMdiff are trademarks of STEMCELL Technologies Canada Inc. iFluor is a trademark of AAT Bioquest, Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.