

Antibodies

Anti-Human CD137 Antibody, Clone 4B4-1, Biotin

Mouse monoclonal IgG1 antibody against human, monkey CD137, biotin-conjugated



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog #100-1354

100 µg

1 mg/mL

Product Description

This mouse monoclonal antibody (clone 4B4-1, biotin-conjugated) reacts with human and monkey CD137, a type I membrane protein and a member of the tumor necrosis factor (TNF) receptor superfamily. Following activation, it is expressed on multiple immune cells, such as natural killer (NK), T, and dendritic cells. CD137 plays a role in cytokine activation, preventing activation-induced cell death, and it also plays a role in promoting the activity of cytotoxic T cells. As an immune stimulator, CD137 can influence the tumor microenvironment by increasing the cytotoxicity of T and NK cells as well as their ability to infiltrate tumors. For these reasons, CD137 represents a strong candidate for cancer immunotherapy. CD137 deficiency in humans has been hypothesized to inhibit differentiation and function of T cells; CD137-deficient mice have been reported to have impaired T cell cytotoxicity, proliferation, and survival.

Target Antigen Name:	CD137
Alternative Names:	4-1BB, CDw137, ILA, TNFRSF9
Gene ID:	3604
Species Reactivity:	Human, Monkey
Host Species:	Mouse
Clonality:	Monoclonal
Clone:	4B4-1
Isotype:	IgG1, kappa
Immunogen:	Fusion protein of recombinant human 4-1BB ectodomain
Conjugate:	Biotin

Applications

Verified:	FC
Reported:	ELISA

Abbreviations: CellSep: Cell separation; ChIP: Chromatin immunoprecipitation; FA: Functional assay; FACS: Fluorescence-activated cell sorting; FC: Flow cytometry; ICC: Immunocytochemistry; IF: Immunofluorescence microscopy; IHC: Immunohistochemistry; IHC-P: Immunohistochemistry (paraffin-embedded); IP: Immunoprecipitation; RIA: Radioimmunoassay; WB: Western blotting

Properties

Formulation:	50 mM sodium phosphate, 100 mM potassium chloride, 150 mM sodium chloride, pH 7.5, containing 5% glycerol, 0.2% (w/v) bovine serum albumin, and 0.04% sodium azide
Purification:	The antibody was purified by affinity chromatography.
Stability and Storage:	Product stable at 2 - 8°C when stored undiluted. Do not freeze. Stable until expiry date (EXP) on label.
Directions for Use:	For flow cytometry, the suggested use of this antibody is 10 µg/mL per 5 x 10 ⁵ cells in 80 µL. It is recommended that the antibody be titrated for optimal performance for each application.

Related Products

For a complete list of antibodies, including other conjugates, sizes, and clones, as well as related products available from STEMCELL Technologies, visit www.stemcell.com/antibodies, or contact us at techsupport@stemcell.com.

References

1. Segal NH et al. (2017) Results from an integrated safety analysis of urelumab, an agonist anti-CD137 monoclonal antibody. *Clinical Cancer Research* 23(8): 1929–36.
2. Somekh I et al. (2019) CD137 deficiency causes immune dysregulation with predisposition to lymphomagenesis. *Blood* 134(18): 1510–6.
3. Seidel J et al. (2021) Role of ADAM10 and ADAM17 in regulating CD137 function. *International Journal of Molecular Sciences* 22(5): 2730.
4. Wu M et al. (2019) Induction of CD137 expression by viral genes reduces T cell costimulation. *Journal of Cellular Physiology* 234(11): 21076–88.
5. Yi L. et al. (2022) CD137 agonists targeting CD137-mediated negative regulation show enhanced antitumor efficacy in lung cancer. *Frontiers in Immunology* 13. (ELISA)
6. Yi L. et al. (2014) Human and mouse CD137 have predominantly different binding CRDs to their respective ligands. *PLOS One* 9(1): e86337. (ELISA)

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.