

**Anti-Human CD137 Antibody,  
Clone 4B4-1, FITC**



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Antibodies

Mouse monoclonal IgG1 antibody against human, monkey CD137, FITC-conjugated

Catalog #100-1355

120 Tests 1.67 µl/test

## Product Description

This mouse monoclonal antibody (clone 4B4-1, FITC-conjugated) reacts with human and monkey CD137, a type I membrane protein and a member of the tumor necrosis factor (TNF) receptor superfamily. Following activation, it is expressed on multiple immune cells, such as natural killer (NK), T, and dendritic cells. CD137 plays a role in cytokine activation, preventing activation-induced cell death, and it also plays a role in promoting the activity of cytotoxic T cells. As an immune stimulator, CD137 can influence the tumor microenvironment by increasing the cytotoxicity of T and NK cells as well as their ability to infiltrate tumors. For these reasons, CD137 represents a strong candidate for cancer immunotherapy. CD137 deficiency in humans has been hypothesized to inhibit differentiation and function of T cells; CD137-deficient mice have been reported to have impaired T cell cytotoxicity, proliferation, and survival.

Target Antigen Name:	CD137
Alternative Names:	4-1BB, CDw137, ILA, TNFRSF9
Gene ID:	3604
Species Reactivity:	Human, Monkey
Host Species:	Mouse
Clonality:	Monoclonal
Clone:	4B4-1
Isotype:	IgG1, kappa
Immunogen:	Fusion protein of recombinant human 4-1BB ectodomain
Conjugate:	FITC (Fluorescein isothiocyanate)

## Applications

Verified:	FC
Reported:	FC

Abbreviations: CellSep: Cell separation; ChIP: Chromatin immunoprecipitation; FA: Functional assay; FACS: Fluorescence-activated cell sorting; FC: Flow cytometry; ICC: Immunocytochemistry; IF: Immunofluorescence microscopy; IHC: Immunohistochemistry; IHC-P: Immunohistochemistry (paraffin-embedded); IP: Immunoprecipitation; RIA: Radioimmunoassay; WB: Western blotting

## Properties

Formulation:	50 mM sodium phosphate, 100 mM potassium chloride, 150 mM sodium chloride, pH 7.5, containing 5% glycerol, 0.2% (w/v) bovine serum albumin, and 0.04% sodium azide
Purification:	The antibody was purified by affinity chromatography.
Stability and Storage:	Product stable at 2 - 8°C when stored undiluted. Do not freeze. Stable until expiry date (EXP) on label.
Directions for Use:	For flow cytometry, the suggested use of this antibody is 1.67 µL per 5 x 10 <sup>5</sup> cells in 80 µL. It is recommended that the antibody be titrated for optimal performance for each application.

## Related Products

For a complete list of antibodies, including other conjugates, sizes, and clones, as well as related products available from STEMCELL Technologies, visit [www.stemcell.com/antibodies](http://www.stemcell.com/antibodies), or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## References

1. Segal NH et al. (2017) Results from an integrated safety analysis of urelumab, an agonist anti-CD137 monoclonal antibody. *Clinical Cancer Research* 23(8): 1929–36.
2. Somekh I et al. (2019) CD137 deficiency causes immune dysregulation with predisposition to lymphomagenesis. *Blood* 134(18): 1510–6.
3. Seidel J et al. (2021) Role of ADAM10 and ADAM17 in regulating CD137 function. *International Journal of Molecular Sciences* 22(5): 2730.
4. Wu M et al. (2019) Induction of CD137 expression by viral genes reduces T cell costimulation. *Journal of Cellular Physiology* 234(11): 21076–88.
5. Wehler T. C. et al. (2008) Rapid identification and sorting of viable virus-reactive CD4+ and CD8+ T cells based on antigen-triggered CD137 expression. *Journal of Immunological Methods* 339(1): 23–37. (FC)

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.