

17 beta-Estradiol

Steroid hormone

Catalog #100-1647

5 g

Product Description

17 beta-Estradiol is the most predominant circulating estrogen secreted by the premenopausal ovary. 17 beta-Estradiol is converted from testosterone primarily in the ovarian granulosa cells and in small amounts in the adrenal gland and brain (Azcoitia et al.; Tsuchiya et al.). 17 beta-Estradiol impacts physiological processes such as reproduction, bone remodeling, and the anti-inflammatory responses by binding to estrogen receptors and regulating gene expression. It has also been implicated to have a role in several diseases including cancer and neurodegenerative disorders (Deroo & Korach).

Alternative Names: 17 β -Oestradiol, β -Estradiol, E2, Estradiol

CAS Number: 50-28-2


Chemical Formula: C₁₈H₂₄O₂

Molecular Weight: 272.4 g/mol

Purity: \geq 98%

Chemical Name: estra-1,3,5(10)-triene-3,17 β -diol

Structure:



Properties

Product Format:	A crystalline solid
Stability and Storage:	Product stable at -20°C as supplied. As a precaution, STEMCELL recommends storing all small molecules away from direct light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
Preparation:	<ul style="list-style-type: none"> • DMSO \leq 70 mM • Absolute ethanol \leq 9.2 mM <p>For example, to prepare a 10 mM stock solution in DMSO, resuspend 10 mg in 3.67 mL of DMSO.</p> <p>Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported; however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.</p> <p>Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO or absolute ethanol concentration above 0.1% due to potential cell toxicity.</p>

Published Applications

IMMUNOLOGY

- Induces immune cell recruitment and a pro-inflammatory response in an influenza mouse model, resulting in improved outcomes (Davis et al.).
- Promotes recruitment of M1-polarized macrophages to inflammation sites in a rat model in a dose-dependent manner (Kou et al.).

CANCER RESEARCH

- Induces epithelial-to-mesenchymal transition activation in human glioblastoma multiforme cells (Hernández-Vega et al.).
- Induces human breast cancer cell cytoskeleton remodeling and enhances migration and invasiveness (Zheng et al.).
- Regulates cell proliferation, migration, and self-renewal capacity of human breast cancer stem cells in a dose-dependent manner (Guo et al.).

References

- Azcoitia I et al. (2011) Estradiol synthesis within the human brain. *Neurosci* 191: 139–47.
- Davis SM et al. (2017) Estradiol and progesterone influence on influenza infection and immune response in a mouse model. *Am J Reprod Immunol* 78(4): e12695.
- Deroo BJ & Korach KS. (2006) Estrogen receptors and human disease. *J Clin Invest* 116(3): 561–70.
- Guo L et al. (2018) 17 β -estradiol regulates the malignancy of cancer stem-like cells derived from the MCF7 cell line partially through Sox2. *Oncol Lett* 15(3): 3790–5.
- Hernández-Vega AM et al. (2020) Estradiol induces epithelial to mesenchymal transition of human glioblastoma cells. *Cells* 9(9): 1930.
- Kou XX et al. (2015) Estradiol promotes M1-like macrophage activation through cadherin-11 to aggravate temporomandibular joint inflammation in rats. *J Immunol* 194(6): 2810–8.
- Tsuchiya Y et al. (2005) Cytochrome P450-mediated metabolism of estrogens and its regulation in human. *Cancer Lett* 227(2): 115–24.
- Zheng S et al. (2011) 17 β -Estradiol enhances breast cancer cell motility and invasion via extra-nuclear activation of actin-binding protein ezrin. *PLoS One* 6(7): e22439.

Related Products

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

Warning

This product is hazardous. Please refer to the Safety Data Sheet (SDS).

17 beta-Estradiol

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.