

UNC0224

Histone modifier; Inhibits euchromatic histone methyltransferase 1 (EHMT1) and 2 (EHMT2)

Catalog #100-1653

5 mg

Product Description

UNC02244 is a potent inhibitor of the euchromatic histone methyltransferase 1 (EHMT1; IC_{50} = 39 nM) and EHMT2 (IC_{50} = 15 nM). EHMT1 and EHMT2 can mono-, di-, and tri-methylate the lysine 9 residue of histone 3 (H3K9) to impact gene regulation (Leenders et al.; Liu et al.). EHMT1 and EHMT2 can heterodimerize to maintain H3K9 methylation in embryonic stem cells (Shinkai & Tachibana; Zhang et al.). Overexpression of EHMT1 and EHMT2 in several cancer types is associated with metastasis, stemness, and therapy resistance (Nachiyappan et al.).

Alternative Names: UNC 0224

CAS Number: 1197196-48-7


Chemical Formula: $C_{26}H_{43}N_7O_2$

Molecular Weight: 485.7 g/mol

Purity: $\geq 98\%$

Chemical Name: 7-[3-(Dimethylamino)propoxy]-2-(hexahydro-4-methyl-1H-1,4-diazepin-1-yl)-6-methoxy-N-(1-methyl-4-piperidiny)-4-quinazolinamine

Structure:



Properties

Product Format:	A white to off-white powder
Stability and Storage:	Product stable at -20°C as supplied. As a precaution, STEMCELL recommends storing all small molecules away from direct light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
Preparation:	<ul style="list-style-type: none">• DMSO \leq 100 mM For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 206 μ L of DMSO. Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported; however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application. Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

Published Applications

REPROGRAMMING

- Reprograms human somatic cells to pluripotent stem cells in combination with other small molecules (Guan et al.).

References

- Guan J et al. (2022) Chemical reprogramming of human somatic cells to pluripotent stem cells. *Nature* 605(7909): 325–31.
- Leenders R et al. (2019) Novel SAR for quinazoline inhibitors of EHMT1 and EHMT2. *Bioorg Med Chem Lett* 29(17): 2516–24.
- Liu F et al. (2009) Discovery of a 2,4-diamino-7-aminoalkoxyquinazoline as a potent and selective inhibitor of histone lysine methyltransferase G9a. *J Med Chem* 52(24): 7950–3.
- Nachiyappan A et al. (2022) EHMT1/EHMT2 in EMT, cancer stemness and drug resistance: emerging evidence and mechanisms. *FEBS J* 289(5): 1329–51.
- Shinkai Y & Tachibana M. (2011) H3K9 methyltransferase G9a and the related molecule GLP. *Genes Dev* 25(8): 781–8.
- Zhang T et al. (2016) G9a/GLP complex maintains imprinted DNA methylation in embryonic stem cells. *Cell Rep* 15(1): 77–85.

Related Products

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.