

ErythroClear[™] Magnet

Magnet for column-free immunomagnetic separation

Catalog #01737

1 Unit

Product Description

ErythroClear[™] Magnet is an anodized aluminum body housing rare earth magnetic elements. ErythroClear[™] Magnet is compatible with procedures using select EasySep[™] reagents when processing up to 16 samples at a time. It is designed to hold up to 16 standard microcentrifuge tubes (e.g. Catalog #38089).

Properties

Stability and Storage: Store at 15 - 25°C. Shelf Life: Not applicable.

Directions for Use

For more information, refer to the appropriate EasySep[™] Product Information Sheet or contact us at techsupport@stemcell.com.

CARE OF MAGNET: After use, wipe the magnet clean with distilled water and then with alcohol to dry the magnet. To disinfect the magnet, wipe with 10% diluted household bleach and then wipe off with a water-dampened cloth to avoid corrosion of metal parts. Do not submerge the magnet in liquid.

Related Products

For a complete list of related products available from STEMCELL Technologies, visit www.stemcell.com, or contact us at techsupport@stemcell.com.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO WWW.STEMCELL.COM/COMPLIANCE.

Copyright© 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, EasySep, and ErythroClear are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.