

Human Recombinant R-spondin-3, ACF

R-spondin-3

Catalog #100-2065

100 µg

Product Description

R-Spondin-3 is a member of thrombospondin type 1 repeat (TSR-1) superfamily that is involved in the canonical Wnt/ β -catenin signaling pathway (de Lau et al.). R-spondin proteins are characterized by two furin-like repeats at the amino terminus and thrombospondin domain located near the carboxyl terminus (de Lau et al.). R-spondin-3 expression is associated with ovarian cancer (Gu et al.), prostate cancer (Mesci et al.), and differentiation of intestinal epithelial cells in diabetes mellitus (Shan et al.). In a transgenic mouse model, the expression of R-Spondin-3 induces the expansion of Lgr5⁺ stem cells, Paneth cells, and Lgr4⁺ cells, promoting the intestinal stem cell compartment (Hilkens et al.). This product is animal component-free (ACF).

Product Information

Alternative Names:	Cristin 1, hRspo3, PWTSR, R-Spondin 3, Roof plate-specific spondin-3, RSPO3, RSpondin 3, Thrombospondin type-1 domain-containing protein 2, THSD2
Accession Number:	Q9BXY4
Predicted Molecular Mass:	17 kDa
Species:	Human
Product Formulation:	Lyophilized from a solution containing acetonitrile and trifluoroacetic acid.
Source:	E. coli
Purity:	≥ 98% by SDS-PAGE

Specifications

Activity:	The EC ₅₀ is approximately 5 ng/mL (~0.3 nM), as determined by the WNT-responsive luciferase reporter assay in HEK239T cells.
Endotoxin Level:	Measured by kinetic Limulus amebocyte lysate (LAL) analysis and is ≤ 0.1 EU/μg protein.

Preparation and Storage

Stability and Storage:	Store at -20 to -80°C. Stable as supplied for 12 months from date of receipt.
Preparation:	<p>Centrifuge vial before opening. Reconstitute the product in 10 mM hydrochloric acid to at least 0.1 mg/mL by pipetting the solution down the sides of the vial. Do not vortex.</p> <p>OPTIONAL: After reconstitution, if product will not be used immediately, dilute with concentrated bovine serum albumin (BSA) to a final BSA concentration of 0.1 - 1%. The effect of storage of stock solution on product performance should be tested for each application. As a general guide, do not store at -20 to -80°C for more than 12 months. Avoid repeated freeze-thaw cycles.</p>

Data




Figure 1. Biological Activity and Molecular Mass of Human Recombinant R-spondin-3, ACF

(A) The biological activity of Human Recombinant R-Spondin-3, ACF was tested by its ability to induce Wnt- β catenin signaling in HEK239T cells using a luciferase reporter assay. Firefly luciferase activity was normalized to control Renilla luciferase activity. The EC₅₀ is defined as the effective concentration of the growth factor at which Wnt- β catenin response is at 50% of maximum. The EC₅₀ in the above example is 0.3 nM (5 ng/mL). (B) 7 μ g of Human Recombinant R-spondin-3, ACF was resolved with SDS-PAGE under reducing (+) and non-reducing (-) conditions and visualized by Coomassie Blue staining. Human Recombinant R-spondin-3, ACF has a predicted molecular mass of 17 kDa.

Related Products

For a complete list of cytokines or peptide pools, as well as related products available from STEMCELL Technologies, visit www.stemcell.com/cytokines or contact us at techsupport@stemcell.com.

References

- Gu H et al. (2020) RSPO3 is a marker candidate for predicting tumor aggressiveness in ovarian cancer. *Ann Transl Med* 8(21): 1351.
- Hilkens J et al. (2017) RSPO3 expands intestinal stem cell and niche compartments and drives tumorigenesis. *Gut* 66(6): 1095–105.
- de Lau WBM et al. (2012) The R-spondin protein family. *Genome Biol* 13(242).
- Mesci A et al. (2019) RSPO3 is a prognostic biomarker and mediator of invasiveness in prostate cancer. *J Transl Med* 17(1): 1–11.
- Shan TD et al. (2021) RSPO3 regulates the abnormal differentiation of small intestinal epithelial cells in diabetic state. *Stem Cell Res Ther* 12(1): 1–11.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.