

EasySep™ Buffer

Cell separation buffer

Catalog #20144

1 L

Product Description

EasySep™ Buffer is recommended for use in most manual EasySep™ cell separation protocols.

Properties

Stability and Storage:

Store at 2 - 8°C. Stable until the expiry date (EXP) on label.

Components may be shipped at room temperature (15 - 25°C) but should be stored as indicated above.

Contains:

- Dulbecco's phosphate-buffered saline (PBS)
- Fetal bovine serum (2%)
- EDTA 1 mM in PBS

Directions for Use

EasySep™ Buffer has been optimized for use with EasySep™ cell separation kits. Refer to the Product Information Sheet (PIS) for the applicable cell separation kit or contact us for more information.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO WWW.STEMCELL.COM/COMPLIANCE.

Copyright © 2025 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and EasySep are trademarks of STEMCELL Technologies Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.