

Anti-Human Beta-Tubulin III Antibody, Clone 2G10-TB3



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Antibodies

Mouse monoclonal IgG2a antibody
against human, mouse, rat beta-
tubulin III, unconjugated

Catalog #60092
#60092.1

100 µg 0.5 mg/mL
25 µg 0.5 mg/mL

Product Description

The 2G10-TB3 antibody reacts with beta-tubulin III, an ~50 - 55 kDa structural protein that is a constituent of tubulin. Tubulin is the major component of microtubules within the cytoskeleton and is assembled from heterodimers of alpha- and beta-tubulin subunits. The beta III isoform of tubulin, also known as neuron-specific class III beta-tubulin, is expressed primarily in neurons and is widely used as a marker to distinguish neurons from other cell types. Beta-tubulin III contributes to microtubule formation in neuronal cell bodies and axons, a function involving GTP binding, and plays roles in axonal transport, neuronal cell proliferation, and differentiation. It is highly expressed in several types of cancer and is a predictive and prognostic marker for various tumors, for example, being found in neoplastic but not in normal glial cells. The 2G10-TB3 antibody is expected to recognize all mammalian homologs of beta-tubulin III and the epitope has reportedly been mapped to the extreme C-terminal amino acids, EAQGPK.

Target Antigen Name:	Beta-Tubulin III
Alternative Names:	Class 3 beta-tubulin, class III beta-tubulin, MC1R, neuron-specific class 3 beta-tubulin, neuron-specific class III beta-tubulin, TUBB 3, TUBB3, tubulin beta 3, tubulin beta 4, tubulin beta III
Gene ID:	10381
Species Reactivity:	Human, Mouse, Rat, Other Mammals
Host Species:	Mouse
Clonality:	Monoclonal
Clone:	2G10-TB3
Isotype:	IgG2a, kappa
Immunogen:	Synthetic peptide corresponding to amino acids 436 - 450 of beta-tubulin III conjugated to keyhole limpet hemocyanin
Conjugate:	Unconjugated

Applications


Verified:	ICC, IF, WB
Reported:	FC, ICC, IF, IHC, IP, WB
Special Applications:	This antibody clone has been verified for labeling neural stem and progenitor cells grown with STEMdiff™ Neural Induction Medium (Catalog #05835), STEMdiff™ Neural Progenitor Medium (Catalog #05833), NeuroCult™ NS-A Proliferation Kit (Human; Catalog #05751), NeuroCult™ Proliferation Kit (Mouse & Rat; Catalog #05702), and NeuroCult™ SM1 Neuronal Culture Kit (Catalog #05712).

Abbreviations: CellSep: Cell separation; ChIP: Chromatin immunoprecipitation; FA: Functional assay; FACS: Fluorescence-activated cell sorting; FC: Flow cytometry; ICC: Immunocytochemistry; IF: Immunofluorescence microscopy; IHC: Immunohistochemistry; IP: Immunoprecipitation; RIA: Radioimmunoassay; WB: Western blotting

Properties

Formulation:	Aqueous buffer containing 0.09% sodium azide, may contain carrier protein/stabilizer
Purification:	The antibody was purified by column chromatography.
Stability and Storage:	Product stable at 2 - 8°C when stored undiluted. Do not freeze. For product expiry date, please contact techsupport@stemcell.com .
Directions for Use:	For immunocytochemistry the suggested use of this antibody is ≤ 5 µg/mL. It is recommended that the antibody be titrated for optimal performance for each application.

Data



(A) E18 cortical rat neurons were cultured using NeuroCult™ SM1 Neuronal Culture Kit on poly-lysine-coated glass coverslips, then fixed and labeled with Anti-Human Beta-Tubulin III Antibody, Clone 2G10-TB3, followed by donkey anti-mouse IgG, Alexa Fluor® 488. Nuclei were counterstained with DAPI. Inset shows cells labeled with a mouse IgG2a, kappa isotype control antibody followed by donkey anti-mouse IgG, Alexa Fluor® 488, and counterstained with DAPI.

(B) Western blot analysis of denatured/reduced cell lysates with Anti-Beta Tubulin III Antibody, Clone 2G10-TB3. Lane 1, adult rat brain cortical cells; lane 2 (negative control), mouse E13.5 neural progenitor cells cultured with NeuroCult™ Proliferation Kit (Mouse & Rat).

Related Products

For a complete list of antibodies, including other conjugates, sizes and clones, as well as related products available from STEMCELL Technologies, please visit our website at www.stemcell.com/antibodies or contact us at techsupport@stemcell.com.

References

1. Nihongaki Y et al. (2017) CRISPR–Cas9-based photoactivatable transcription systems to induce neuronal differentiation. *Nat Methods* 14(10): 963–6. (IF)
2. Malchenko S et al. (2014) Onset of rosette formation during spontaneous neural differentiation of hESC and hiPSC colonies. *Gene* 534(2): 400–7. (FC, ICC, IF)
3. Karki R et al. (2013) β III-Tubulin: biomarker of taxane resistance or drug target? *Expert Opin Ther Targets* 17(4): 461–72.
4. Locher H et al. (2013) Class III β -tubulin, a novel biomarker in the human melanocyte lineage. *Differentiation* 85(4-5): 173–81. (ICC, IF, IHC, WB)
5. Mobarakeh ZT et al. (2012) Human endometrial stem cells as a new source for programming to neural cells. *Cell Biol Int Rep* 19(1): e00015. (ICC, IF)
6. Mariani M et al. (2011) Class III β -tubulin (TUBB3): more than a biomarker in solid tumors? *Curr Mol Med* 11(9): 726–31.
7. Yu H et al. (2011) Lentiviral gene transfer into the dorsal root ganglion of adult rats. *Mol Pain* 7: 63. (ICC, IF, IHC, WB)
8. Koh Y et al. (2009) Class III beta-tubulin, but not ERCC1, is a strong predictive and prognostic marker in locally advanced head and neck squamous cell carcinoma. *Ann Oncol* 20(8): 1414–9. (IHC)
9. Dráberová E et al. (1998) Expression of class III beta-tubulin in normal and neoplastic human tissues. *Histochem Cell Biol* 109(3): 231–9. (ICC, IF, IHC, WB)
10. Ludueña RF. (1998) Multiple forms of tubulin: different gene products and covalent modifications. *Int Rev Cytol* 178: 207–75.
11. Lee MK et al. (1990) Posttranslational modification of class III beta-tubulin. *Proc Natl Acad Sci USA* 87(18): 7195–9. (Immunoaffinity chromatography, WB)

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, STEMdiff, and NeuroCult are trademarks of STEMCELL Technologies Canada Inc. Matrigel is a trademark of Corning® Incorporated. All other trademarks are the property of their respective holders. Alexa Fluor is a registered trademark of Life Technologies Corporation. Antibodies conjugated to Alexa Fluor® are licensed for internal research use only and sale is expressly conditioned on the buyer not using the antibody for manufacturing, performing a service or medical test, or otherwise generating revenue. For use other than research, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@lifetech.com. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.