

# Small Molecules

## Dorsomorphin

BMP and AMPK pathway inhibitor;  
Inhibits ALK2, ALK3, ALK6, and AMPK

Catalog # 72102

10 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

Dorsomorphin inhibits the bone morphogenetic protein (BMP) pathway by targeting the type I BMP receptors activin receptor-like kinase (ALK) 2, ALK3, and ALK6. It is also a potent inhibitor of AMP-activated protein kinase (AMPK,  $K_i = 109$  nM) but does not significantly inhibit structurally related kinases such as ZAPK, SYK, PKC $\theta$ , PKA, or JAK3 (Bain et al.; Yu et al.).

Molecular Name:	Dorsomorphin
Alternative Names:	Compound C
CAS Number:	866405-64-3
Chemical Formula:	C <sub>24</sub> H <sub>25</sub> N <sub>5</sub> O
Molecular Weight:	399.5 g/mol
Purity:	≥ 98%
Chemical Name:	6-[4-[2-(1-piperidinyl)ethoxy]phenyl]-3-(4-pyridinyl)-pyrazolo[1,5-a]pyrimidine
Structure:	



## Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at -20°C as supplied. Protect from prolonged exposure to light. Stable as supplied for 12 months from date of receipt.
Solubility:	· DMSO ≤ 12 mM (with heat applied) · Absolute ethanol ≤ 350 μM · Dimethylformamide (DMF) ≤ 6.2 mM For example, to prepare a 300 μM stock solution in absolute ethanol, resuspend 1 mg in 8.3 mL of absolute ethanol.  Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in absolute ethanol at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.  Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final ethanol concentration above 0.1% due to potential cell toxicity.

## Published Applications

### DIFFERENTIATION

- Promotes differentiation of neural progenitor cells from human pluripotent stem cells (Morizane et al.; Zhou et al.).
- Promotes differentiation of cardiomyocytes from mouse and human pluripotent stem cells (Hao et al.; Kattman et al.).
- Promotes differentiation of adipocytes and suppresses osteogenic differentiation of osteoblasts from human mesenchymal cells (Kim et al.).

## References

- Bain J et al. (2007) The selectivity of protein kinase inhibitors: a further update. *Biochem J* 408(3): 297–315.
- Hao J et al. (2008) Dorsomorphin, a selective small molecule inhibitor of BMP signaling, promotes cardiomyogenesis in embryonic stem cells. *PLoS One* 3(8): e2904.
- Kattman SJ et al. (2011) Stage-specific optimization of activin/nodal and BMP signaling promotes cardiac differentiation of mouse and human pluripotent stem cell lines. *Cell Stem Cell* 8(2): 228–40.
- Kim E-K et al. (2012) Human mesenchymal stem cell differentiation to the osteogenic or adipogenic lineage is regulated by AMP-activated protein kinase. *J Cell Physiol* 227(4): 1680–7.
- Morizane A et al. (2011) Small-molecule inhibitors of bone morphogenic protein and activin/nodal signals promote highly efficient neural induction from human pluripotent stem cells. *J Neurosci Res* 89(2): 117–26.
- Yu PB et al. (2008) Dorsomorphin inhibits BMP signals required for embryogenesis and iron metabolism. *Nat Chem Biol* 4(1): 33–41.
- Zhou J et al. (2010) High-efficiency induction of neural conversion in human ESCs and human induced pluripotent stem cells with a single chemical inhibitor of transforming growth factor beta superfamily receptors. *Stem Cells* 28(10): 1741–50.

## Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit [www.stemcell.com/smallmolecules](http://www.stemcell.com/smallmolecules) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

**This product is hazardous. Please refer to the Safety Data Sheet (SDS).**

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.