

Anti-R-Phycoerythrin (R-PE) Antibody, Clone ID3



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog #01402 100 µL 1 mg/mL

Product Description

R-phycoerythrin (R-PE) can be coupled to antibodies, peptides, or ligands to allow their detection by fluorescent microscopy or flow cytometry. The ID3 antibody reacts with R-PE alone or with PE-tandem dyes, such as PE-Cyanine 7, without quenching the fluorescent signal. When combined into tetrameric antibody complexes (TACs), the ID3 antibody can be used to isolate highly pure populations of PE-labeled cells.

This antibody clone is used in EasySep™ PE Positive Selection Kits (Catalog #18551 [Human], #18554 [Mouse], and #18557) and in EasySep™ Release PE Positive Selection Kits (Catalog #17654 [Human] and #17656 [Mouse]) to isolate cells of interest at high purity from single-cell suspensions derived from human, mouse, or other species. Additionally, it can be used with EasySep™ Human “Do-It-Yourself” Selection Kit (Catalog #18099) to build custom TACs for positive selection from human single-cell suspensions.

Target Antigen Name:	R-Phycoerythrin (R-PE)
Species Reactivity:	Not applicable
Host Species:	Mouse (BALB/c)
Clonality:	Monoclonal
Clone:	ID3
Isotype:	IgG1, kappa
Immunogen:	Purified R-phycoerythrin protein
Conjugate:	Unconjugated

Applications

Verified:	CellSep
Reported:	ELISA, FC

Properties

Formulation:	Phosphate-buffered saline
Purification:	The antibody was purified by affinity chromatography.
Stability and Storage:	Product stable at 2 - 8°C until expiry date (EXP) as indicated on label. Do not freeze. Addition of 0.1% sodium azide (final) is recommended once the vial has been opened.
Directions for Use:	Centrifuge vial briefly before use to ensure recovery of contents. The antibody should be titrated for optimal performance for each application.

Related Products

For a complete list of antibodies, including other conjugates, sizes, and clones, as well as related products available from STEMCELL Technologies, visit www.stemcell.com/antibodies or contact us at techsupport@stemcell.com.

References

1. Lansdorp PM et al. (1991) Single laser three color immunofluorescence staining procedures based on energy transfer between phycoerythrin and cyanine 5. *Cytometry* 12(8): 723–30. (ELISA, FC)
2. Wognum AW et al. (1987) Use of tetrameric antibody complexes to stain cells for flow cytometry. *Cytometry* 8(4): 366–71.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and EasySep are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.