

L-Glutamine

Cell culture supplement

Catalog # 07100 100 mL



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

L-Glutamine (200 mM) is recommended for use as a cell culture supplement and is required for the preparation of MesenCult™-XF Medium (Catalog #05420/05429) and MesenCult™-ACF Medium (Catalog #05440/05449).

Properties

Storage: Store at -20°C.
Shelf Life: Stable until expiry date (EXP) on label.

Handling / Directions For Use

1. Thaw L-Glutamine in a 37°C water bath. Mix well.
NOTE: If not used immediately, aliquot and store at -20°C.
2. Add 5 mL of L-Glutamine to 500 mL of culture medium to achieve a final concentration of 2 mM.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, MesenCult, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.