



## RoboSep™ Tip Head Polishing Compound

Catalog #20119



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

The RoboSep™ Tip Head Polishing Compound is a cleaning/polishing solution designed to remove residue buildup that may accumulate over time, thus improving the tip contact surface condition for pickup/stripping of disposable RoboSep™ Tips.

## Contains

- 70% Isopropanol
- Tween 20

## Stability and Storage

Product stable at 15 - 25°C.

## Directions For Use

Apply 1 - 2 drops of RoboSep™ Tip Head Polishing Compound to a laboratory wipe. Rub onto the RoboSep™ tip head. With a fresh laboratory wipe, wipe off and dry the tip head. Repeat weekly.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and RoboSep are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.