

40 µm Cell Strainer

40 µm nylon cell strainer

Catalog # 27305

50 Strainers/Case



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

These sterile 40 µm nylon cell strainers are useful for obtaining single-cell suspensions during enzymatic digestion of mammary tissue and organoids. They may also be used for the dissociation of cells from other primary tissue, and in the preparation of single-cell suspensions for flow cytometry.

Properties

Storage: Store at 15 - 25°C.

Shelf Life: Not applicable.

Handling / Directions For Use

NOTE: For single use only.

Unwrap cell strainer unit, fit into a 50 mL conical tube (e.g. Catalog #38010), and filter cells through the strainer.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.