

# Blunt-End Needles, 16 Gauge



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## 16 gauge needles with luer lock

|                 |           |
|-----------------|-----------|
| Catalog # 28110 | 100/Pack  |
| 28120           | 2000/Pack |

## Product Description

Sterile 16 gauge needles, 1.5 inches long with a metal luer lock hub are ideal for safe and accurate handling of methylcellulose-based medium (e.g. MethoCult™). The blunt end reduces the risk of needle-stick injuries. Product is manufactured from non-pyrogenic materials.

## Properties

**Storage:** Store at 15 - 25°C.  
**Shelf Life:** Stable until expiry date (EXP) on label.

## Handling / Directions For Use

Use a 3 cc Syringe (Catalog #28230) and a Blunt-End Needle to dispense MethoCult™ methylcellulose-based medium. This is for accurate dispensing of viscous methylcellulose medium and to prevent needle-stick injuries.

Use Blunt-End Needles with syringes purchased from STEMCELL Technologies to ensure that they are compatible. STEMCELL Technologies cannot guarantee that Blunt-End Needles will be compatible with syringes from other suppliers.

Blunt-End Needles and 3 cc Syringes are single-use only. Dispose according to standard procedures for biohazardous materials.

For complete instructions, refer to one of the following Technical Manuals, available at [www.stemcell.com](http://www.stemcell.com) or contact us to request a copy:

- Human Colony-Forming Unit Assays Using MethoCult™ (Document #28404)
- Mouse Colony-Forming Unit Assays Using MethoCult™ (Document #28405)

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2019 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and MethoCult are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.