

Iscove's Modified Dulbecco's Medium



Iscove's Modified Dulbecco's Medium (IMDM) with 25 mM HEPES

Catalog # 36150 500 mL

Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Iscove's Modified Dulbecco's Medium (IMDM) is a synthetic medium that is recommended for a wide variety of cell culture applications. Selection of a suitable nutrient medium is dependent on the type of cell, conditions of culture, and degree of chemical definition required for the cell culture application.

See Formulation on page 2 for a full list of components.

Properties

Storage: Store at 2 - 8°C.
Shelf Life: Stable for 12 months from date of manufacture (MFG) on label.
Contains: IMDM with 25 mM HEPES

Handling / Directions For Use

This product contains the labile amino acid L-glutamine which has a half-life of approximately one month when stored at 2 - 8°C. For certain cell culture applications, medium stored for greater than two months following the date of manufacture should be supplemented with additional L-glutamine. Add 5 mL of 200 mM L-Glutamine (Catalog #07100) to 500 mL to achieve a final concentration of 2 mM.

NOTE: This product does not contain antibiotics. If desired, add penicillin and streptomycin and use medium within 1 week.

Formulation

INORGANIC SALTS	g/L
Calcium Chloride•2H ₂ O	0.219
Magnesium Sulfate	0.09767
Potassium Chloride	0.33
Potassium Nitrate	0.000076
Sodium Chloride	4.505
Sodium Phosphate Monobasic	0.109
Sodium Selenite	0.000017
AMINO ACIDS	g/L
L-Alanine	0.025
L-Arginine•HCl	0.084
L-Asparagine•H ₂ O	0.0284
L-Aspartic Acid	0.03
L-Cystine•2HCl	0.09124
L-Glutamic Acid	0.075
L-Glutamine	0.584
Glycine	0.03
L-Histidine•HCl•H ₂ O	0.042
L-Isoleucine	0.105
L-Leucine	0.105
L-Lysine•HCl	0.146
L-Methionine	0.03
L-Phenylalanine	0.066
L-Proline	0.04
L-Serine	0.042
L-Threonine	0.095
L-Tryptophan	0.016
L-Tyrosine 2Na•2H ₂ O	0.10379
L-Valine	0.094
VITAMINS	g/L
D-Biotin	0.000013
Choline Chloride	0.004
Folic Acid	0.004
myo-Inositol	0.0072
Niacinamide	0.004
D-Pantothenic Acid (hemicalcium)	0.004
Pyridoxine•HCl	0.004
Riboflavin	0.0004
Thiamine•HCl	0.004
Vitamin B-12	0.000013
OTHER	g/L
D-Glucose	4.5
HEPES	5.958
Phenol Red•Na	0.016
Pyruvic Acid•Na	0.11
Sodium Bicarbonate	3.024

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2015 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Inc. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.