

Small Molecules

Gö6983

PKC inhibitor

Catalog # 72462

1 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Gö6983 inhibits several isoforms of protein kinase C (PKC; IC_{50} = 7, 7, 6, 10, 60, and 20,000 nM for PKC α , PKC β , PKC γ , PKC δ , PKC ζ , and PKC μ , respectively; Gschwendt et al.).

Molecular Name:	Gö6983
Alternative Names:	Gö 6983
CAS Number:	133053-19-7
Chemical Formula:	C ₂₆ H ₂₆ N ₄ O ₃
Molecular Weight:	442.5 g/mol
Purity:	≥ 98%
Chemical Name:	3-[1-[3-(dimethylamino)propyl]-5-methoxy-1H-indol-3-yl]-4-(1H-indol-3-yl)-1H-pyrrole-2,5-dione
Structure:	



Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at -20°C as supplied. Protect from prolonged exposure to light. Stable as supplied for 12 months from date of receipt.
Solubility:	· DMSO ≤ 20 mM For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 226 μ L of fresh DMSO.

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

Published Applications

MAINTENANCE AND SELF-RENEWAL

- Inhibits differentiation and maintains pluripotency in mouse embryonic stem cells (Dutta et al.).
- Enhances human ground state pluripotent stem cell viability and growth (Gafni et al.).
- Inhibits proliferation of human primary fetal bone cells (Krattinger et al.).
- Inhibits the formatin of proplatelets from megakaryocytes derived from adult mouse bone marrow (Williams et al.).

REPROGRAMMING

- Enhances reprogramming of mouse embryonic fibroblasts to induced pluripotent stem cells (Dutta et al.).
- Direct lineage reprogramming of fibroblasts to mature neurons, in combination with CHIR99021 (Catalog #72052), RepSox (Catalog #73792), Forskolin (Catalog #72112), SP600125 (Catalog #72642), Valproic Acid (Catalog #72292), and Y-27632 (Catalog #72302) (Hu et al.).

References

- Dutta D et al. (2011) Self-renewal versus lineage commitment of embryonic stem cells: protein kinase C signaling shifts the balance. *Stem Cells* 29(4): 618–28.
- Gafni O et al. (2013) Derivation of novel human ground state naive pluripotent stem cells. *Nature* 504(7479): 282–6.
- Gschwendt M et al. (1996) Inhibition of protein kinase C μ by various inhibitors. Inhibition from protein kinase c isoenzymes. *FEBS Lett* 392(2): 77–80.
- Hu W et al. (2015) Direct conversion of normal and Alzheimer's Disease human fibroblasts into neuronal cells by small molecules. *Cell Stem Cell* 17(2): 204–12.
- Krattinger N et al. (2011) Regulation of proliferation and differentiation of human fetal bone cells. *Eur Cell Mater* 21: 46–58.
- Williams CM et al. (2014) PKC α negatively regulates in vitro proplatelet formation and in vivo platelet production in mice. *Platelets* 25(1): 62–8.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

This product is hazardous. Please refer to the Safety Data Sheet (SDS).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.