

Small Molecules

PS-48

PI3K/AKT pathway activator; Activates PDK1

Catalog # 72832
72834

10 mg
50 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

PS-48 is a phosphoinositide-dependent protein kinase-1 (PDK1) activator. It acts by binding to the HM/PIF binding pocket of PDK1 rather than the ATP-binding site ($K_d = 10.3 \mu\text{M}$; Hindie et al.; Stroba et al.). PDK1 is a master kinase upstream of the AGC kinases (AKT, S6K, RSK, SGK, and PKC) in response to increases in the second messenger PtdIns-(3,4,5)-P₃ (Bayascas).

Molecular Name:	PS-48
Alternative Names:	Not applicable
CAS Number:	1180676-32-7
Chemical Formula:	C ₁₇ H ₁₅ ClO ₂
Molecular Weight:	286.8 g/mol
Purity:	≥ 98%
Chemical Name:	(Z)-5-(4-chlorophenyl)-3-phenylpent-2-enoic acid
Structure:	



Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at -20°C as supplied. Protect from prolonged exposure to light. For product expiry date, please contact techsupport@stemcell.com.
Solubility:	<ul style="list-style-type: none">· Absolute ethanol ≤ 100 mM· DMSO ≤ 100 mM For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 349 μL of fresh DMSO.

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

Published Applications

REPROGRAMMING

· In combination with A-83-01 and PD0325901, enhances reprogramming of neonatal human epidermal keratinocytes transduced with OCT4 and KLF4 (Zhu et al.).

References

Bayascas JR. (2010) PDK1: the major transducer of PI 3-kinase actions. *Curr Top Microbiol Immunol* 346: 9–29.

Hindie V et al. (2009) Structure and allosteric effects of low-molecular-weight activators on the protein kinase PDK1. *Nat Chem Biol* 5(10): 758–64.

Stroba A et al. (2009) 3,5-Diphenylpent-2-enoic acids as allosteric activators of the protein kinase PDK1: structure-activity relationships and thermodynamic characterization of binding as paradigms for PIF-binding pocket-targeting compounds. *J Med Chem* 52(15): 4683–93.

Zhu S et al. (2010) Reprogramming of human primary somatic cells by OCT4 and chemical compounds. *Cell Stem Cell* 7(6): 651–5.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, please visit our website at www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2015 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design and Scientists Helping Scientists are trademarks of STEMCELL Technologies Inc. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.