

Cytokines

Rat Recombinant GM-CSF



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Granulocyte-macrophage colony-stimulating factor

Catalog #	78018.1	20 µg
	78018	100 µg
	78018.2	1000 µg

Product Description

Granulocyte-macrophage colony-stimulating factor (GM-CSF) promotes the proliferation and differentiation of hematopoietic progenitor cells and the generation of neutrophils, eosinophils, and macrophages. In synergy with other cytokines such as stem cell factor, IL-3, erythropoietin, and thrombopoietin, it also stimulates erythroid and megakaryocyte progenitor cells (Barreda et al.). GM-CSF is produced by multiple cell types, including stromal cells, Paneth cells, macrophages, dendritic cells (DCs), endothelial cells, smooth muscle cells, fibroblasts, chondrocytes, and Th1 and Th17 T cells cells (Francisco-Cruz et al.). The receptor for GM-CSF (GM-CSFR) is composed of two subunits: the cytokine-specific α subunit (GMR α ; CD116) and the common subunit β c (CD131) shared with IL-3 and IL-5 receptors (Broughton et al.). GM-CSFR is expressed on hematopoietic cells, including progenitor cells and immune cells, as well as non-hematopoietic cells. GM-CSF is able to stimulate the development of DCs that ingest, process, and present antigens to the immune system (Francisco-Cruz et al.). Recombinant rat GM-CSF is reactive with mouse cells (Oaks et al.; Vandenabeele et al.).

Product Information

Alternative Names:	Colony-stimulating factor 2, CSF-2, MGI-1GM, Pluripoietic-alpha
Accession Number:	P48750
Amino Acid Sequence:	MAPTRSPNPV TRPWKHVDAI KEALLLNDM RALENEKNED VDIISNEFSI QRPTCVQTRL KLYKQGLRGN LTKLNGALTM IASHYQTNCP PTPETDCEIE VTFEDFIKN LKGFLFDIPF DCWKPVQK
Predicted Molecular Mass:	14.7 kDa
Species:	Rat
Cross Reactivity:	Mouse
Formulation:	Lyophilized from a sterile-filtered aqueous solution containing sodium bicarbonate, pH 8.5.
Source:	E. coli

Specifications


Activity:	The specific activity is $\geq 5 \times 10^7$ units/mg ($EC_{50} \leq 20$ pg/mL) as determined by a cell proliferation assay using FDC-P1 cells.
Purity:	$\geq 95\%$
Endotoxin Level:	Measured by kinetic Limulus amoebocyte lysate (LAL) analysis and is ≤ 1 EU/ μ g protein.

Preparation and Storage

Storage:	Store at -20°C to -80°C.
Stability:	Stable as supplied for 12 months from date of receipt.
Preparation:	Centrifuge vial before opening. Reconstitute the product in sterile water to at least 0.1 mg/mL by pipetting the solution down the sides of the vial. Do not vortex.

OPTIONAL: After reconstitution, if product will not be used immediately, dilute with concentrated bovine serum albumin (BSA) to a final BSA concentration of 0.1%. The effect of storage of stock solution on product performance should be tested for each application. As a general guide, do not store at 2 - 8°C for more than 1 month or at -80°C for more than 3 months. Avoid repeated freeze-thaw cycles.

Data



(A) The biological activity of Rat Recombinant GM-CSF was tested by its ability to promote the proliferation of FDC-P1 cells. Cell proliferation was measured after 91 hours of culture using a fluorometric assay method. The EC₅₀ is defined as the effective concentration of the growth factor at which cell proliferation is at 50% of maximum. The EC₅₀ in the above example is 0.9 - 1.3 pg/mL. (B) 1 µg of Rat Recombinant GM-CSF was resolved with SDS-PAGE under reducing (+) and non-reducing (-) conditions and visualized by Coomassie Blue staining. Rat Recombinant GM-CSF has a predicted molecular mass of 14.7 kDa.

Related Products

For a complete list of cytokines, as well as related products available from STEMCELL Technologies, visit www.stemcell.com/cytokines or contact us at techsupport@stemcell.com.

References

- Barreda DR et al. (2004) Regulation of myeloid development and function by colony stimulating factors. *Dev Comp Immunol* 28(5): 509–54.
- Broughton SE et al. (2012) The GM-CSF/IL-3/IL-5 cytokine receptor family: from ligand recognition to initiation of signaling. *Immunol Rev* 250(1): 277–302.
- Francisco-Cruz A et al. (2014) Granulocyte-macrophage colony-stimulating factor: not just another haematopoietic growth factor. *Med Oncol* 31(1): 774.
- Oaks MK et al. (1995) Polymerase chain reaction cloning and expression of the rat granulocyte-macrophage colony-stimulating factor. *J Interferon Cytokine Res* 15(12): 1095–102.
- Vandenabeele P et al. (1990) Response of murine cell lines to an IL-1/IL-2-induced factor in a rat/mouse T hybridoma (PC60): differential induction of cytokines by human IL-1 alpha and IL-1 beta and partial amino acid sequence of rat GM-CSF. *Lymphokine Res* 9(3): 381–9.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.