

**Anti-Human MUC1 (CD227)
Antibody, Clone 214D4**

Antibodies

Mouse monoclonal IgG1 antibody
against human MUC1 (CD227),
unconjugated

Catalog #60137

200 µg 1 mg/mL



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

The 214D4 antibody reacts with MUC1 (also known as CD227), a large cell surface glycoprotein expressed by most glandular and ductal epithelial cells as well as a variety of hematopoietic cells. A characteristic feature of the MUC1 glycoprotein is a core protein domain composed of a variable number of tandem repeats and multiple oligosaccharide side chains. Because the extracellular portion of MUC1 can extend beyond most cell surface proteins, it is thought to play a role in cell-cell and cell-substrate adhesion. The protein is highly expressed by the majority of human adenocarcinomas and is associated with a poor prognosis. In the mammary gland, MUC1 is localized on the apical plasma membrane of luminal epithelial cells. The 214D4 antibody clone is particularly useful since it enriches for luminal-restricted but not bipotent epithelial progenitor cells from normal human mammary tissue.

Target Antigen Name:	MUC1 (CD227)
Alternative Names:	CD227, EMA, Episialin, Epithelial membrane antigen, HMFG antigen, MAM6, Mucin 1, PEM, Polymorphic epithelial mucin
Gene ID:	4582
Species Reactivity:	Human
Host Species:	Mouse
Clonality:	Monoclonal
Clone:	214D4
Isotype:	IgG1
Immunogen:	Recombinant protein comprising the mucin-like repeat domain of human MUC1 (CD227) fused to beta-galactosidase
Conjugate:	Unconjugated

Applications


Verified:	FC, IHC
Reported:	FC, IHC

Abbreviations: CellSep: Cell separation; ChIP: Chromatin immunoprecipitation; FA: Functional assay; FACS: Fluorescence-activated cell sorting; FC: Flow cytometry; ICC: Immunocytochemistry; IF: Immunofluorescence microscopy; IHC: Immunohistochemistry; IP: Immunoprecipitation; RIA: Radioimmunoassay; WB: Western blotting

Properties

Formulation:	Tris-glycine, pH 7.4, containing sodium chloride and 0.05% sodium azide before the addition of glycerol to 30%
Purification:	The antibody was purified by protein G chromatography.
Stability and Storage:	Product stable at -20°C when stored undiluted. For product expiry date, please contact techsupport@stemcell.com.
Directions for Use:	The suggested use of this antibody is: FC, 1 - 3 µg/mL per 1 x 10 ⁶ cells in 100 µL; IHC, 1 - 3 µg/mL. It is recommended that the antibody be titrated for optimal performance for each application.

Data



Flow cytometry analysis of human MCF7 cells labeled with Anti-Human MUC1 (CD227) Antibody, Clone 214D4, followed by a rat anti-mouse IgG1 antibody, FITC (filled histogram) or a mouse IgG1, kappa isotype control antibody, followed by a rat anti-mouse IgG1 antibody, FITC (solid line histogram).

Related Products

For a complete list of antibodies, including other conjugates, sizes and clones, as well as related products available from STEMCELL Technologies, please visit our website at www.stemcell.com/antibodies or contact us at techsupport@stemcell.com.

References

1. Silva AL et al. (2014) MUC1 expression in Fallopian tubes of women with hydrosalpinx. *Eur J Obstet Gynecol Reprod Biol* 180: 106–10. (IHC)
2. Martignani E et al. (2010) Human milk protein production in xenografts of genetically engineered bovine mammary epithelial stem cells. *PLoS One* 5(10): e13372. (FC)
3. Eirew P et al. (2008) A method for quantifying normal human mammary epithelial stem cells with in vivo regenerative ability. *Nat Med* 14(12): 1384–9. (IHC)
4. O'Connor JC et al. (2005) MUC1 expression in human prostate cancer cell lines and primary tumors. *Prostate Cancer Prostatic Dis* 8(1): 36–44.
5. Jarrard JA et al. (1998) MUC1 is a novel marker for the type II pneumocyte lineage during lung carcinogenesis. *Cancer Res* 58(23): 5582–9.

Please refer to the Safety Data Sheet (SDS) for hazard information.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.