

## Small Molecules

### Anacardic Acid

Epigenetic modifier; Inhibits histone acetyltransferase (HAT)

Catalog # 73192  
73194

1 mg  
5 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

Anacardic Acid is a 6-alkyl salicylic acid that inhibits the histone acetyltransferase (HAT) activity of the transcription co-activators p300 and p300/CREB-binding protein-associated factor (pCAF; IC<sub>50</sub> values of 8.5 and 5 μM, respectively; Balasubramanyam et al.). Anacardic Acid inhibits HAT-dependent transcription and protein SUMOylation (Cui et al.; Fukuda et al.). In addition, Anacardic Acid is an activator of Aurora kinase A-mediated phosphorylation of Histone H3 (Kishore et al.). At higher concentrations (25 μM), Anacardic Acid suppresses NF-κB activation and inhibits IκB-α phosphorylation (Sung et al.).

Molecular Name:	Anacardic Acid
Alternative Names:	Not applicable
CAS Number:	16611-84-0
Chemical Formula:	C <sub>22</sub> H <sub>36</sub> O <sub>3</sub>
Molecular Weight:	348.5 g/mol
Purity:	≥ 98%
Chemical Name:	2-hydroxy-6-pentadecylbenzoic acid
Structure:	



## Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at -20°C as supplied. Protect from prolonged exposure to light. For long-term storage store with a desiccant. Stable as supplied for 12 months from date of receipt.
Solubility:	· DMSO ≤ 25 mM · Absolute ethanol ≤ 25 mM For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 287 μL of DMSO.

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

## Published Applications

### DIFFERENTIATION

- Enhances cardiomyocyte differentiation from mouse embryonic stem cells (Re et al.).

### IMMUNOLOGY

- Induces macrophage activation (Gnanaprakasam et al.).
- Exhibits inhibitory and bactericidal activities against methicillin-resistant *Staphylococcus aureus* (Muroi & Kubo).

### CANCER RESEARCH

- Potentiates the apoptosis induced by cytokine and chemotherapeutic agents in cancer cells (Sung et al.).
- Sensitizes tumor cells to ionizing radiation in vitro (Sun et al.).

## References

- Balasubramanyam K et al. (2003) Small molecule modulators of histone acetyltransferase p300. *J Biol Chem* 278(21): 19134–40.
- Cui L et al. (2008) Histone acetyltransferase inhibitor anacardic acid causes changes in global gene expression during in vitro *Plasmodium falciparum* development. *Eukaryot Cell* 7(7): 1200–10.
- Fukuda I et al. (2009) Ginkgolic acid inhibits protein SUMOylation by blocking formation of the E1-SUMO intermediate. *Chem Biol* 16(2): 133–40.
- Gnanaprakasam JN et al. (2015) The anacardic 6-pentadecyl salicylic acid induces macrophage activation via the phosphorylation of ERK1/2, JNK, P38 kinases and NF- $\kappa$ B. *Int Immunopharmacol* 29(2): 808–17.
- Kishore AH et al. (2008) Specific small-molecule activator of Aurora kinase A induces autophosphorylation in a cell-free system. *J Med Chem* 51(4): 792–7.
- Muroi H & Kubo I. (1996) Antibacterial activity of anacardic acid and totarol, alone and in combination with methicillin, against methicillin-resistant *Staphylococcus aureus*. *J Appl Bacteriol* 80(4): 387–94.
- Re A et al. (2015) Anacardic acid and thyroid hormone enhance cardiomyocytes production from undifferentiated mouse ES cells along functionally distinct pathways. *Endocrine* 53(3): 681–8.
- Sun Y et al. (2006) Inhibition of histone acetyltransferase activity by anacardic acid sensitizes tumor cells to ionizing radiation. *FEBS Lett* 580(18): 4353–6.
- Sung B et al. (2008) Anacardic acid (6-nonadecyl salicylic acid), an inhibitor of histone acetyltransferase, suppresses expression of nuclear factor- $\kappa$ B-regulated gene products involved in cell survival, proliferation, invasion, and inflammation through inhibition of the inhibitory subunit of nuclear factor- $\kappa$ B $\alpha$  kinase, leading to potentiation of apoptosis. *Blood* 111(10): 4880–91.

## Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit [www.stemcell.com/smallmolecules](http://www.stemcell.com/smallmolecules) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

**This product is hazardous. Please refer to the Safety Data Sheet (SDS).**

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.