

## Small Molecules

### Butein

NF- $\kappa$ B pathway inhibitor; Stabilizes I $\kappa$ B

Catalog # 73462  
73464

1 mg  
10 mg



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

Butein is a plant polyphenol, tetrahydroxychalcone that inhibits nuclear factor (NF)- $\kappa$ B (Yang et al; Pandey et al. 2007). Butein has been shown to prevent phosphorylation and degradation of tumor necrosis factor (TNF)-dependent I $\kappa$ B $\alpha$ , an inhibitory subunit of NF- $\kappa$ B (IC<sub>50</sub> = 38  $\mu$ M; Orlikova et al.). Butein also inhibits 5-lipoxygenase (IC<sub>50</sub> = 0.01  $\mu$ M; Sogawa et al.), the enoyl-acyl-carrier protein reductase of *Plasmodium falciparum* (Ki = 2.97  $\mu$ M; Sharma et al.), angiotensin-converting enzyme (IC<sub>50</sub> = 0.73 mM; Bonesi et al.), and SRC kinase (Pandey et al. 2009).

Molecular Name:	Butein
Alternative Names:	2',3,4,4'-Tetrahydroxychalcone
CAS Number:	487-52-5
Chemical Formula:	C <sub>15</sub> H <sub>12</sub> O <sub>5</sub>
Molecular Weight:	272.3 g/mol
Purity:	≥ 95%
Chemical Name:	(2E)-1-(2,4-dihydroxyphenyl)-3-(3,4-dihydroxyphenyl)-2-propen-1-one
Structure:	



## Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at -20°C as supplied. Protect from prolonged exposure to light. For long-term storage store with a desiccant. Stable as supplied for 12 months from date of receipt.
Solubility:	· DMSO ≤ 70 mM · Ethanol ≤ 70 mM For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 367 $\mu$ L of DMSO.

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

## Published Applications

### METABOLISM

- Inhibits iron-induced lipid peroxidation in rat brain homogenate in a concentration-dependent manner (Cheng et al.).

### CANCER RESEARCH

- Inhibits TNF- $\alpha$ -induced invasion of human lung adenocarcinoma H1299 cells (Pandey et al. 2007; Gupta et al.).
- Inhibits growth in human hepatoma cancer cell lines HepG2 and Hep3B, by inducing G2/M phase arrest (Moon et al.; Gupta et al.).
- Inhibits testosterone-induced proliferation in breast cancer cells (Wang et al.).

### IMMUNOLOGY

- Exhibits anti-inflammatory properties in a mouse macrophage cell line by inhibiting lipopolysaccharide-induced expression of inducible nitric oxide synthase (Lee et al.).

## References

- Bonesi M et al. (2010) The synthesis and angiotensin converting enzyme (ACE) inhibitory activity of chalcones and their pyrazole derivatives. *Bioorg Med Chem Lett* 20(6): 1990–3.
- Cheng ZJ et al. (1998) Antioxidant properties of butein isolated from *Dalbergia odorifera*. *Biochim Biophys Acta* 1392(2-3): 291–9.
- Gupta SC et al. (2010) Regulation of survival, proliferation, invasion, angiogenesis, and metastasis of tumor cells through modulation of inflammatory pathways by nutraceuticals. *Cancer Metastasis Rev* 29(3): 405–34.
- Lee SH et al. (2004) Inhibition of lipopolysaccharide-induced expression of inducible nitric oxide synthase by butein in RAW 264.7 cells. *Biochem Biophys Res Commun* 323(1): 125–32.
- Orlikova B et al. (2012) Natural chalcones as dual inhibitors of HDACs and NF- $\kappa$ B. *Oncol Rep* 28(3): 797–805.
- Moon D-O et al. (2010) Butein induces G(2)/M phase arrest and apoptosis in human hepatoma cancer cells through ROS generation. *Cancer Lett* 288(2): 204–13.
- Pandey MK et al. (2007) Butein, a tetrahydroxychalcone, inhibits nuclear factor (NF)- $\kappa$ B and NF- $\kappa$ B-regulated gene expression through direct inhibition of I $\kappa$ B kinase beta on cysteine 179 residue. *J Biol Chem* 282(24): 17340–50.
- Pandey MK et al. (2009) Butein suppresses constitutive and inducible signal transducer and activator of transcription (STAT) 3 activation and STAT3-regulated gene products through the induction of a protein tyrosine phosphatase SHP-1. *Mol Pharmacol* 75(3): 525–33.
- Sharma SK et al. (2007) Green tea catechins potentiate triclosan binding to enoyl-ACP reductase from *Plasmodium falciparum* (PfENR). *J Med Chem* 50(4): 765–75.
- Sogawa S et al. (1993) 3,4-Dihydroxychalcones as potent 5-lipoxygenase and cyclooxygenase inhibitors. *J Med Chem* 36(24): 3904–9.
- Wang Y et al. (2005) The plant polyphenol butein inhibits testosterone-induced proliferation in breast cancer cells expressing aromatase. *Life Sci* 77(1): 39–51.
- Yang EB et al. (1998) Butein, a specific protein tyrosine kinase inhibitor. *Biochem Biophys Res Commun* 245(2): 435–8.

## Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit [www.stemcell.com/smallmolecules](http://www.stemcell.com/smallmolecules) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.